

SARS-COV-2 ANTIGEN TEST KIT

(colloidal gold method)


BioTeke

USER INSTRUCTION



1. Read this instruction guide carefully.
2. Prepare a watch (or a clock/timer), tissues and either hand sanitizer or soap and warm water.
3. Check the test kit contents. Make sure that nothing is damaged or broken.

- For anterior nasal swabs.
- For in-vitro self-testing use only.
- Please read the instructions carefully before you begin testing!



1

Wash your hands thoroughly for at least 20 seconds before the test.




2

Put the Tube into the kit box holder before proceeding to the tube. Add all of Sample Diluent vertically to the Sample extraction tube.

3


NOTE: Please blow your nose before collection.

Remove the swab from its wrapper and take out the swab by holding the handle. Being careful not to touch the fabric tip of the swab with your hands.



4

Gently insert the swab into your nostril less than one inch (about 2.5cm). Slowly rub the swab against all of the inside walls of your nostril. Make at least 5 big circles. Do not just spin the swab. Repeat this step in your other nostril using the same swab.



NOTE: With children, the maximum depth of insertion into the nostril may be less than 3/4 of an inch.

5

Insert the swab into the tube. Touch the bottom of the tube with the swab tip, and stir at least 15 times. Squeeze the swab in the tube through the outer wall of the tube by finger 5 times.



6

Screw on the purple cap and unscrew the top white cap.



7

Open the pouch and take out the Test Card. Place it on a flat, dry and clean surface. Turn the tube integrated dropper cap upside down and slowly squeeze 4 drops onto the sample well of the Test Card.



15-20 minutes


8 Results Interpretation



NOTE: The test results should not be read after 20 minutes.

COVID-19 Detected (Positive)

A positive test result indicates that antigens from COVID-19 were detected, and the patient is very likely to be infected with the virus and presumed to be contagious. Test results should always be considered in the context of clinical observations and epidemiological data in making a final diagnosis and patient management decisions. You should self-isolate at home and avoid contact with others as per CDC recommendations to avoid spreading the virus to others.



COVID-19 Not Detected (Negative)

A negative test result indicates that antigens from COVID-19 were not detected from the specimen. However, a negative result does not rule out COVID-19 and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of an individual's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19 and confirmed with a molecular assay, if necessary, for patient management.



Invalid barcode or absence of a purple-colored line next to "C".

Re-test with a COVID-19 test may be needed. An invalid test result indicates that your test has experienced an error and is unable to interpret the result of the test. You will need to re-test with a new test or consult a healthcare professional. If you still have symptoms, you should self-isolate at home and avoid contact with others prior to the retest.



9

All used test components should be disposed of in your household waste. After completing all steps, wash hands or use hand sanitizer.

