Dräger Antigen Test SARS-CoV-2

Instructions for use

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de	Gebrauchsanweisung.....	3
en	Instructions for use	11
fr	Notice d'utilisation.....	19
es	Instrucciones de uso.....	28
it	Istruzioni per l'uso.....	36
nl	Gebruiksaanwijzing	44

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13 Literaturhinweise

- (1) Lindner, A. K. et al. 2021 Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab. *Eur. Respir. J.* 57: 2003961 doi:10.1183/13993003.03961-2020.
- (2) Sungnak, W. et al. 2020, SARS-CoV-2 entry factors are highly expressed in nasal epithelial cells together with innate immune genes. *Nature Medicine* 26, 681–687. doi.org/10.1038/s41591-020-0868-6
- (3) Robert Koch-Institut, 2021. COVID-19: Entlassungskriterien aus der Isolierung. Robert Koch-Institut. doi.org/10.25646/6957.6
- (4) Wölfel, R. et al., 2020. Virological assessment of hospitalized patients with COVID-2019. *Nature* 581, 465–469. doi.org/10.1038/s41586-020-2196-x
- (5) van Kampen, J.J.A. et al., 2021. Duration and key determinants of infectious virus shedding in hospitalized patients with coronavirus disease-2019 (COVID-19). *Nat Commun* 12, 267. doi.org/10.1038/s41467-020-20568-4
- (6) Pollock, N.R. et al., 2021. Correlation of SARS-CoV-2 nucleocapsid antigen and RNA concentrations in nasopharyngeal samples from children and adults using an ultrasensitive and quantitative antigen assay. *J Clin Microbiol* 59(4) doi.org/10.1128/JCM.03077-20
- (7) Osmanodja, B. et al., 2021. Accuracy of a Novel SARS-CoV-2 Antigen-Detecting Rapid Diagnostic Test from Standardized Self-Collected Anterior Nasal Swabs. *J Clin Med* 10(10), 2099 <https://doi.org/10.3390/jcm10102099>

⚠ Read the instructions for use carefully before carrying out the test. Carry out each step in the order described. In addition, take a look at the video.

The following QR code will take you to the video:

<http://antigentest.draeger.com/use>

1 Intended use

The Antigen Test SARS-CoV-2 is a test kit consisting of an immunochemical rapid test based on test strips and a sample collector for sampling the anterior nasal cavity. As a representative virus component, the antigen test qualitatively detects the nucleocapsid protein (antigen) of the human SARS-CoV-2 virus from the nasal swab of persons suspected of having a SARS-CoV-2 infection.

The antigen test is intended for professional use by healthcare professionals only.

The antigen test provides a preliminary result (screening test), which needs to be confirmed using a further measurement method. A PCR test can be used for this purpose, for example. The results of the antigen test require a medical evaluation taking into account all relevant clinical information about the patient.

2 Background information

Coronaviruses form a large family of viruses that can cause simple colds, but also serious respiratory diseases. In the case of the human coronavirus 2 (SARS-CoV-2), the resulting "COVID-19" disease is harmless in the vast majority of people, but it can also lead to severe acute respiratory syndrome (SARS).

When someone is infected with a virus, human cells reproduce all virus components. This includes the target antigen of this rapid test, the nucleocapsid protein, which is detectable in the upper respiratory tract during the acute phase of infection. For more information, see Chapter 13 References. ^(1 - 2)

3 Test principle

The antigen test is based on the principle of the sandwich immunoassay. Viral protein (here nucleocapsid protein) present in the sample is bound by both capture antibodies and red labelled detection antibodies (gold conjugate).

The reaction of the antibodies with the viral protein is shown as a pink line in area T (test line) of the viewing window. Depending on the virus concentration of the sample, the intensity of this line may vary.

3.1 Quality control

Each test strip contains an additional function control independent of the nasal swab. The formation of a pink line in area C (control line) of the viewing window indicates that the test membrane has been completely wetted and that the antibody-antigen-reaction has proceeded as expected.

3.2 Target antigen and calibrator

- SARS-CoV-2 nucleocapsid protein

3.3 Effective reagents

- Gold conjugate: Mouse anti SARS-CoV-2 nucleocapsid protein IgG
- Test line: Mouse anti SARS-CoV-2 nucleocapsid protein IgG
- Control line: Goat anti-Chicken IgY

4 Safety-related information

4.1 Basic safety instructions

- Follow the instructions for use carefully. The user must fully understand the instructions and follow them precisely. Use the product only for the purposes specified in the intended use section of this document.
- Do not dispose of the instructions for use. Ensure that they are stored and used appropriately by the product user.
- Only trained or expert persons may use this product.
- Do not use a faulty or incomplete product. Do not modify the product.
- Notify Dräger in the event of any component fault or failure.
- Comply with all local and national rules and regulations associated with this product.

4.2 Product-related safety instructions

- Wear personal protective equipment when handling samples (gloves, spectacles, protective mask, protective clothing).
- Observe legal safety instructions regarding handling of infectious materials.
- Observe decontamination and disposal regulations for infectious materials.
- Do not open the foil pouch until immediately before use.

- Only use the antigen test if the foil pouch is undamaged.
- Do not use the antigen test if the expiration date on the foil pouch has expired. The expiration date is in the format YYYY-MM. Example: 2022-01 means that the antigen test must not be used after the end of January 2022.
- ☒ Test kits must not be reused.
- Treat all samples as if they contained infectious substances.
- Do not smoke, eat or drink during the test procedure.

5 Conventions in this document

The following warnings are used in this document to inform the user of possible dangers. The meanings of the warnings are defined as follows:

Warning sign	Signal word	Classification of the warning
	CAUTION	Indication of a potentially dangerous situation. If this is not avoided, injuries can occur. Can also be used as a warning against improper use.


6 Description

6.1 Product overview

6.1.1 Dräger Antigen Test SARS-CoV-2


Test kit consists of:

- 1 Sample collector
- 2 Test cassette
- 3 Tear-away security seal over viewing window with test strip
- 4 Indicator ring
- 5 Funnel-shaped test opening with ampoule inside (not visible)
- 6 Desiccant bag (not needed for the test)

6.1.2 Viewing window with test strip

- 7 Control line (C)
- 8 Test line (T)

Each packaging unit contains:

- 20 test kits, individually packed in a sealed foil pouch
- 1 instructions for use

7 Use

⚠ CAUTION


Risk of contamination!

Contamination may occur if insufficient protective equipment is worn and the usual hygiene rules are not observed.

- ▶ Wear protective clothing, gloves, spectacles and a protective mask.
- ▶ Do not touch the white end of the sample collector before and after sampling.

⚠ Read the instructions for use carefully before carrying out the test. Carry out each step in the order described. In addition, take a look at the video.

The following QR code will take you to the video:

<http://antigentest.draeger.com/use>

7.1 Sampling

Prerequisites

- The patient must not have applied any substances, such as nasal spray, to the nose at least 15 minutes before sampling.
- Before using the antigen test, check that the test is safe to use on the patient due to the nature of the nose and the inner walls of the nose.
- The ambient temperature and the temperature of the test kit must be +15 to +30 °C.
- Use the test kit immediately after opening the foil pouch.
- If the test kit is dropped, do not use it again.
- The desiccant bag of the test kit must be present.
- The foil pouch must not be damaged.
- The indicator ring must be completely coloured pink before sampling.

Work equipment

- Disposable gloves for handling infectious material (e.g., latex or nitrile gloves)
- Spectacles
- Protective clothing
- Protective mask
- Timer
- Tissue
- Optional: resealable plastic pouch

1. Tear open the foil pouch.
2. Sampling with the antigen test is done by the patient under professional guidance and supervision (**steps 1 to 4**). Instruct the patient to blow their nose vigorously into a tissue.

⚠ Only wipe the nose externally so that sufficient nasal secretion remains in the nostrils (Fig. 1).

ℹ If it is not possible for the patient to take the sample, the healthcare professionals will do so. (steps 1 to 4)


1

48690

3. Pull the sample collector against the resistance and out of the test kit (Fig. 2).


2

48651

4. Insert the white end of the sample collector into a nostril up to the nasal constriction and wipe the inner walls of the nose for 30 seconds in a circular motion using gentle pressure (Fig. 3). Using the same sample collector, repeat the sample collection in the other nostril for 30 seconds.


3


48619

5. The patient hands over the sample collector to the professional user who carries out all further testing steps.
6. Insert the sample collector into the funnel-shaped test opening of the test cassette and press it down firmly against resistance as far as it will go to break the ampoule (Fig 4).

4

5


48652

⚠ Make sure that there is no gap between the sample collector and the test cassette (Fig. 5).

7.2 Sample preparation

1. Shake the test kit well for at least 30 seconds until the indicator ring has visibly discoloured (Fig. 6).


6

48653


2. Then place the test kit on a flat surface for 1 minute (Fig. 7).

7


48655


3. After the wait time has elapsed, completely tear off and remove the security seal on the test kit (**Fig. 8**).

4. In an upright position on a flat surface, forcefully press the test kit down onto the sample collector as far as it will go (**Fig. 9**). The test solution will then flow upwards along the test strip. If the test solution does not flow along the test strip within one minute, tap the test kit firmly on a flat surface.


5. Place the test kit back on a flat surface (**Fig. 10**). As soon as the test solution starts to flow upwards along the test strip, set the timer to 15 minutes and activate.

7.3 Test result

The test result must be read after **15 to 20 minutes at the latest** under adequate lighting.

Positive test result Negative test result Invalid test result

Positive test result

If a pink test line (T) appears, this indicates a positive test result. The intensity of the test line may vary. Thus, even a very weak, partially discoloured or interrupted test line can be considered **positive**.

Negative test result

Negative test results must be read 15 minutes at the earliest and **no later** than 20 minutes after starting the test. If no test line has appeared during this time, while the control line is present, the test is **negative**.

Test is invalid

If neither a pink control line nor a pink test line appears within 20 minutes of starting the test (Chapter 7.2 Sample preparation, Fig. 9), the test is **invalid**.

7.4 Hygienic protective measures (optional)

When using a commercially available, transparent, resealable plastic pouch (not included in the scope of delivery), sample preparation and reading of the test result can be carried out under special hygienic protective precautions, as described in Chapters 7.2 and 7.3.


⚠ CAUTION

- ▶ Put the test kit into the plastic pouch and close the plastic pouch.
- ▶ Make sure that the test kit is placed upright in the closed plastic pouch on a flat surface.

Prerequisites

- Use an unused plastic pouch
- Only use highly transparent, colourless plastic pouches

1. Put the test kit into the plastic pouch and close it.

2. All further steps and the evaluation can be done inside the plastic pouch, as described in Chapters 7.2 and 7.3.

-  After starting the test, place the test kit in the plastic bag upright on a flat surface.
- Do not reuse plastic pouches.
- Provide sufficient lighting when reading the test result.

7.5 Restriction of the procedure

With immunochemical tests, factors may occur that affect the reliability. In very rare cases, this can lead to a positive test result although the patient does not have the SARS-CoV-2 virus (false positive result). For this reason, a more comprehensive confirmation analysis (e.g., PCR) is indispensable to confirm all positive test results.

In rare cases, false negative results may also occur. This means that the antigen test shows a negative result, whereas a reference method (e.g. PCR) gives a positive result. This can especially occur in the early stages of the infection, a few days after contagion or during the decay phase of the infection.

Positive test results do not rule out infection with other pathogens. Negative test results do not rule out a SARS-CoV-2 infection.

8 Performance characteristics

8.1 Diagnostic performance characteristics

The sensitivity and specificity of the Dräger Antigen Test SARS-CoV-2 were determined in two prospective studies during the period from February to March 2021 in Berlin (Germany) and Kiev (Ukraine).

The test results were compared with commercial PCR tests based on a combined oronasopharyngeal swab (Berlin study) or oropharyngeal swab (Kiev study).

The sensitivity of the test was evaluated using swabs from subjects who showed Covid-19 symptoms within 7 days of symptom onset.

Berlin study

Positive PCR samples of the Berlin study were sequenced to identify the viral variant. The Dräger Antigen Test SARS-CoV-2 detected variant B.1.1.7 in addition to the wild type.

Table 1 Berlin study: Results (symptomatic subjects)

	PCR ¹⁾			
	Positive	Negative	Total	
Dräger Antigen Test SARS-CoV-2	Positive	62	1	63
	Negative	7	202	209
	Total	69	203	272

1) Roche Cobas SARS-CoV-2

Table 2 Berlin study: Sensitivity depending on CT value

CT value ¹⁾	Sensitivity	95 % CI ²⁾
All (≤40)	89.9 %	80.5 – 95.0 %
<30	93.9 %	85.4 - 97.6 %
<26	96.8 %	89.0 - 99.1 %

- 1) Target sequence E-Gen
- 2) Confidence interval

A sensitivity of 96.8 % (95 % CI:89.0 – 99.1 %) was achieved with a CT value <26. All samples above CT 26 contained less than 10⁶ RNA copies/mL.

The concentration of 10⁶ RNA copies/mL represents an orientation value for infectivity. ⁽³⁻⁵⁾

Specificity was determined on the basis of all study results that were negative in the PCR tests.

Table 3 Berlin study: Specificity based on symptomatic and asymptomatic subjects

	PCR	
	Positive	Negative ¹⁾
Dräger Antigen Test SARS-CoV-2	Positive	1
	Negative	307
	Total	308
Specificity	99.7 % (95 % CI: 98.2 – 99.9 %)	

1) This data set contains 105 PCR-negative samples from asymptomatic subjects.

The performance evaluation test showed a specificity of 99.7 % (95 % CI: 98.2 – 99.9 %) across all PCR-negative samples.

Kiev study

No type identification of the SARS-CoV-2 viruses was carried out as part of the Kiev study.

Table 4 Kiev study: Results (symptomatic subjects)

	PCR ¹⁾			
	Positive	Negative	Total	
Dräger Antigen Test SARS-CoV-2	Positive	71	1	72
	Negative	12	100	112
	Total	83	101	184

1) Seegene Allplex SARS-CoV-2 Assay

Table 5 Kiev study: Sensitivity depending on CT value

CT value ¹⁾	Sensitivity	95 % CI ²⁾
All (≤40)	85.5 %	75.5 - 91.9 %
<30	90.9 %	81.6 - 96.0 %
<26	95.5 %	86.6 - 98.8 %

- 1) Target sequence E-Gen
- 2) Confidence interval

A sensitivity of 90.9 % (95 % CI:81.6 - 96.0 %) was achieved with a CT value <30 (here, the orientation value for a virus concentration that is no longer infectious).

Table 6 Kiev study: Specificity based on symptomatic and asymptomatic subjects

		PCR Negative ¹⁾
Dräger Antigen Test SARS-CoV-2	Positive	2
	Negative	383
	Total	385
Specificity	99.5 % (95 % CI: 97.9 - 99.9 %)	

1) This data set contains 284 PCR-negative samples from asymptomatic subjects.

Specificity based on all PCR-negative cases was 99.5 % (95.0 % CI 97.9 - 99.9 %).

Study synopsis

The following tables summarize case counts and performance characteristics from both studies. In total, 845 patient samples were evaluated, of which 468 were obtained by healthcare professionals and 377 by self-sampling.

Table 7 Sensitivity synopsis (symptomatic subjects)

CT value ¹⁾	PCR positive	Sensitivity	95 % CI ²⁾
All (≤40)	152	87.5 %	81.3 - 91.9 %
<30	143	92.3 %	86.8 - 95.7 %
<26	129	96.1 %	91.3 - 98.3 %

1) Target sequence E-Gen
2) Confidence interval

A sensitivity of 92.3 % (95 % CI: 86.8 – 95.7 %) was achieved with a CT value <30.

Table 8 Specificity synopsis

		PCR Negative ¹⁾
Dräger Antigen Test SARS-CoV-2	Positive	3
	Negative	690
	Total	693
Specificity	99.6 % (95 % CI: 98.7 - 99.9 %)	

1) This data set contains 389 PCR-negative samples from asymptomatic subjects.

The combined specificity was 99.6 % (95 % CI: 98.7 – 99.9 %).

8.2 Analytical performance characteristics

8.2.1 Detection limit

The detection limits of the Dräger Antigen Test SARS-CoV-2 were determined with dilutions of inactivated virus particles (ATCC 52287), as well as with the target antigen (nucleocapsid protein).

Reference material	Detection limit
Nucleocapsid protein	200 pg/test
ATCC 52287	7 x 10 ³ TCID ₅₀ /mL

8.2.2 Repeatability and reproducibility

Precision was determined by testing 10 replicates of negative, weak and strong positive controls.

Reproducibility was determined by testing negative, weak and strong positive controls in quintuplicate. The influence of users, antigen test batches and light conditions was investigated.

The Dräger Antigen Test SARS CoV-2 showed acceptable repeatability and reproducibility. The negative and positive samples were correctly determined in 100 % of the cases.

8.2.3 Cross-reactivities

The antigen test was examined for possible cross-reactivity to other respiratory pathogens. The organisms specified in the table led to negative test results.

Pathogen	Test concentration
Coxsackie virus B	1 x 10 ⁵ TCID ₅₀ /mL
Human coronavirus 229E	1 x 10 ⁵ TCID ₅₀ /mL
Human coronavirus NL63	1 x 10 ⁵ TCID ₅₀ /mL
Human coronavirus OC43	1 x 10 ⁵ TCID ₅₀ /mL
Human coronavirus HKU1	1 µg nucleocapsid protein/test
Human adenovirus type 1	1 x 10 ⁵ TCID ₅₀ /mL
Human metapneumovirus	1 x 10 ⁵ TCID ₅₀ /mL
Human parainfluenza virus 1	1 x 10 ⁵ TCID ₅₀ /mL
Human parainfluenza virus 2	1 x 10 ⁵ TCID ₅₀ /mL
Human parainfluenza virus 3	1 x 10 ⁵ TCID ₅₀ /mL
Human parainfluenza virus 4	1 x 10 ⁵ TCID ₅₀ /mL
Influenza A	1 x 10 ⁵ TCID ₅₀ /mL
Influenza B	1 x 10 ⁵ TCID ₅₀ /mL
MERS coronavirus	1 x 10 ⁵ TCID ₅₀ /mL
Respiratory syncytial virus	1 x 10 ⁵ TCID ₅₀ /mL
Rhinovirus A	1 x 10 ⁵ TCID ₅₀ /mL
Rhinovirus B	1 x 10 ⁵ TCID ₅₀ /mL
Haemophilus influenzae	1 x 10 ⁶ pfu/mL
Legionella pneumophila	1 x 10 ⁶ pfu/mL
Mycoplasma pneumoniae	1 x 10 ⁶ pfu/mL
Staphylococcus aureus	1 x 10 ⁶ pfu/mL
Streptococcus pneumoniae	1 x 10 ⁶ pfu/mL
Streptococcus pyogenes	1 x 10 ⁶ pfu/mL

8.2.4 Potentially interfering substances

In both negative and positive samples (0.5 - 5 ng nucleocapsid protein/test), no interfering influences by the substances listed in the table were detected.

Substance	Concentration
Mucin	1300 mg/mL
Blood (human)	5 % (v/v)
Phenylephrine (neo-synephrine)	10 % (v/v)
Dexpanthenol in 0.5 % sodium chloride solution	10 % (v/v)
Beclometasone dipropionate (Rhinivict)	10 % (v/v)
Oxymetazoline (Nasivin)	10 % (v/v)
EMSER salt (saline nasal gel)	5 % (v/v)
Zanamivir (Relenza)	5 mg/mL
Oseltamivir (Tamiflu)	10 mg/mL
Ciprofloxacin	30 µM
Doxycycline	70 µM
Ephedrine	100 ng/mL
Amphetamine	100 ng/mL
Cocaine	100 ng/mL
Metamphetamine	100 ng/mL
Tetrahydrocannabinol (THC)	30 ng/mL
Ibuprofen	2.5 mM
Acetylsalicylic acid	3.7 mM

8.2.5 Hook effect

No hook effect was detected in a highly positive sample (8100 ng nucleocapsid protein/test).⁽⁶⁾

9 Alarm – Cause – Remedy

Error	Cause	Remedy
No control lines appear.	Test solution does not run over test strips.	Tap test kit bottom firmly on a flat surface.
	Antigen test does not stand on a flat surface.	
	Antigen test was not started in a vertical position.	
No test lines appear by mistake.	The test strip has no contact with the sample liquid.	Insert the sample collector into the funnel-shaped test opening of the test cassette and press it down firmly against a resistance as far as it will go to break the ampoule.
	Test kit too deep after sampling, pressed directly into end position (as under Fig. 9).	Cancel test, repeat with new test kit.

10 Storage

The antigen test must be stored at temperatures between +4 °C and +25 °C.

11 Disposal

Dispose of all sample and test materials in health care facilities in accordance with local regulations for the disposal of biological hazardous materials.

12 Explanation of type-identifying marking and symbols

Symbol	Explanation
	Follow the instructions for use.
	Only intended for single use
	Temperature range of storage conditions
	This product complies with the requirements of the European Directive 98/79/EC
	Batch
	Order number
	Medical device for in-vitro diagnostics
	Do not use if packaging is damaged!
	Caution
	Expiration date
	Manufacturer
	Contents sufficient for <n> tests

13 References

- (1) Lindner, A. K. et al. 2021 Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab. *Eur. Respir. J.* 57: 2003961 doi:10.1183/13993003.03961-2020.
- (2) Sungnak, W. et al. 2020, SARS-CoV-2 entry factors are highly expressed in nasal epithelial cells together with innate immune genes. *Nature Medicine* 26, 681–687. doi.org/10.1038/s41591-020-0868-6
- (3) Robert Koch-Institut, 2021. COVID-19: Entlassungskriterien aus der Isolierung. Robert Koch-Institut. doi.org/10.25646/6957.6
- (4) Wölfel, R. et al., 2020. Virological assessment of hospitalized patients with COVID-2019. *Nature* 581, 465–469. doi.org/10.1038/s41586-020-2196-x
- (5) van Kampen, J.J.A. et al., 2021. Duration and key determinants of infectious virus shedding in hospitalized patients with coronavirus disease-2019 (COVID-19). *Nat Commun* 12, 267. doi.org/10.1038/s41467-020-20568-4
- (6) Pollock, N.R. et al., 2021. Correlation of SARS-CoV-2 nucleocapsid antigen and RNA concentrations in nasopharyngeal samples from children and adults using an ultrasensitive and quantitative antigen assay. *J Clin Microbiol* 59(4) doi.org/10.1128/JCM.03077-20
- (7) Osmanodja, B. et al., 2021. Accuracy of a Novel SARS-CoV-2 Antigen-Detecting Rapid Diagnostic Test from Standardized Self-Collected Anterior Nasal Swabs. *J Clin Med* 10(10), 2099 <https://doi.org/10.3390/jcm10102099>

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13 Bibliografia

- (1) Lindner, A. K. et al. 2021 Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab. Eur. Respir. J. 57: 2003961 doi:10.1183/13993003.03961-2020.
- (2) Sungnak, W. et al. 2020, SARS-CoV-2 entry factors are highly expressed in nasal epithelial cells together with innate immune genes. Nature Medicine 26, 681–687. doi.org/10.1038/s41591-020-0868-6
- (3) Robert Koch-Institut, 2021. COVID-19: Entlassungskriterien aus der Isolierung. Robert Koch-Institut. doi.org/10.25646/6957.6
- (4) Wölfel, R. et al., 2020. Virological assessment of hospitalized patients with COVID-2019. Nature 581, 465–469. doi.org/10.1038/s41586-020-2196-x
- (5) van Kampen, J.J.A. et al., 2021. Duration and key determinants of infectious virus shedding in hospitalized patients with coronavirus disease-2019 (COVID-19). Nat Commun 12, 267. doi.org/10.1038/s41467-020-20568-4
- (6) Pollock, N.R. et al., 2021. Correlation of SARS-CoV-2 nucleocapsid antigen and RNA concentrations in nasopharyngeal samples from children and adults using an ultrasensitive and quantitative antigen assay. J Clin Microbiol 59(4) doi.org/10.1128/JCM.03077-20
- (7) Osmanodja, B. et al., 2021. Accuracy of a Novel SARS-CoV-2 Antigen-Detecting Rapid Diagnostic Test from Standardized Self-Collected Anterior Nasal Swabs. J Clin Med 10(10), 2099 https://doi.org/10.3390/jcm10102099

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- (1) Lindner, A. K. et al. 2021 Head-to-head comparison of SARS-CoV-2 antigen-detecting rapid test with self-collected anterior nasal swab versus professional-collected nasopharyngeal swab. Eur. Respir. J. 57: 2003961 doi:10.1183/13993003.03961-2020.
- (2) Sungnak, W. et al. 2020, SARS-CoV-2 entry factors are highly expressed in nasal epithelial cells together with innate immune genes. Nature Medicine 26, 681–687. doi.org/10.1038/s41591-020-0868-6
- (3) Robert Koch-Institut, 2021. COVID-19: Entlassungskriterien aus der Isolierung. Robert Koch-Institut. doi.org/10.25646/6957.6
- (4) Wölfel, R. et al., 2020. Virological assessment of hospitalized patients with COVID-2019. Nature 581, 465–469. doi.org/10.1038/s41586-020-2196-x
- (5) van Kampen, J.J.A. et al., 2021. Duration and key determinants of infectious virus shedding in hospitalized patients with coronavirus disease-2019 (COVID-19). Nat Commun 12, 267. doi.org/10.1038/s41467-020-20568-4
- (6) Pollock, N.R. et al., 2021. Correlation of SARS-CoV-2 nucleocapsid antigen and RNA concentrations in nasopharyngeal samples from children and adults using an ultrasensitive and quantitative antigen assay. J Clin Microbiol 59(4) doi.org/10.1128/JCM.03077-20
- (7) Osmanodja, B. et al., 2021. Accuracy of a Novel SARS-CoV-2 Antigen-Detecting Rapid Diagnostic Test from Standardized Self-Collected Anterior Nasal Swabs. J Clin Med 10(10), 2099 <https://doi.org/10.3390/jcm10102099>

 Manufacturer
Dräger Safety AG & Co. KGaA
Revalstraße 1
D-23560 Lübeck
Germany
+49 451 8 82-0

 Directive 98/79/EC

9300366 – 4755.240 me
© **Dräger Safety AG & Co. KGaA**
Edition: 5 – 2021-06 (Edition: 1 – 2020-10)
Subject to alterations
www.draeger.com