IMMUNOBIO

SARS-CoV-2 Antigen Rapid Test

BfArM Test-ID: AT228/20 (naso-/oropharyngeal) | AT150/21 (anterior nasal) | AT149/21 (Speichel)

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EC Declaration of Conformity

Manufacturer:

Name: Hangzhou Immuno Biotech Co.,Ltd.
Address: Room401, 4F, Building 2, No.28, No.3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China.
Tel: 0571-85368996
Email: info@jgbiotech.com

Whose Authorized Representative:

Name: Lotus NL B.V.
Address: Koningin Julianaplein 10,1e Verd, 2595AA, The Hague, Netherlands.
E-mail: peter@lotusnl.com

We, the manufacturer, here with declare that the product(s)

Product Name	SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag)	Specification	40 T/kit, 20 T/kit, 10 T/kit, 1 T/kit.
Intended Use	The SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is a rapid chromatographic immunoassay for the qualitative detection of novel coronavirus SARS-CoV-2 in human throat secretions, nasal secretions and saliva specimen.		
Classification	Others		

Conformity Assessment Route : IVDD98/79/EC Annex III.

Applicable Standards:

- | | | |
|----------------------------|----------------------------|------------------------|
| <i>ISO 13485:2016</i> | <i>EN ISO 18113-3:2011</i> | <i>EN 13612:2002</i> |
| <i>ISO 14971:2019</i> | <i>EN 13641:2002</i> | <i>ISO 23640:2015</i> |
| <i>EN ISO 18113-1:2011</i> | <i>ISO 15223-1:2016</i> | <i>EN 62366-1:2015</i> |
| <i>EN ISO 18113-2:2011</i> | | |

We, the manufacturer, here declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We agree to develop, implement and maintain a documented post-production monitoring process.

Name of General Manager	Xiaoguang Li
Signature	
Date	Nov 30 2020
Place	Hangzhou, China.
Seal (Manufacturer)	

This is to certify that the Quality Management System of

Hangzhou Immuno Biotech Co., Ltd.

Unified Social Credit Code: 91330101092042370C

Operation Address: Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

Registered Address: Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

applicable to

R&D, Production and Sales of Coronavirus COVID-19 IgG/IgM Antibody Rapid Test, SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag), SARS-CoV-2 Neutralizing Antibody Rapid Test (COVID-19 Ab)(Export to EU)

has been assessed and registered by NQA against the provisions of

ISO 13485: 2016

This registration is subject to the company maintaining a quality management system, to the above standard, which will be monitored by NQA.

Certified Clients shall accept regular surveillance assessments, the validity of certificates shall be maintained for the positive result of audit.

The information of this certificate can be checked on CNCA's website (www.cnca.gov.cn)

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Managing Director

Certificate Number **47938**

Date: 27 July 2020
Reissue Date: 13 January 2021
Valid Until: 27 July 2023

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Hangzhou Immuno Biotech Co., Ltd.

Adr.: Room 401, 4F, Building 2, No. 28, No. 3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China
Tel.: 0571-85368996
E-Mail: info@jgbiotech.com

EU-Bevollmächtigter:

Lotus NL B.V.

Koningin Julianaplein 10,1e Verd,2595AA, The Hague, Niederlande.
E-Mail: Peter@lotusnl.com

IMMUNOBIO


SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19) is a rapid chromatographic immunoassay for the qualitative detection of novel coronavirus SARS-CoV-2 in human nasal cavity.

Hangzhou Immuno Biotech Co.,Ltd.
Room401, 4F, Building 2, No.28, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

European Authorized Representative: Lotus NL B. V.
Koninqln Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands.
E-mail: Peter@lotusnl.com

CONTENT OF THE KIT

PRECAUTIONS FOR USE


- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS



- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
- *NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).

SPECIMEN COLLECTION

TEST PROCEDURE

IMMUNOBIO


SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is a rapid test for the qualitative detection of antigens to the novel coronavirus SARS-CoV-2 in sputum.

Hangzhou Immuno Biotech Co., Ltd.
Room 401, 4F, Building 2, No.26, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

European Authorized Representative: Lotus NL B. V.
Koningin Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands.
E-mail: Peter@lotusnl.com

CONTENT OF THE KIT

PRECAUTIONS FOR USE

- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS

- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
- *NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).

SPECIMEN COLLECTION


Put the specimen collection tube(a) on the rack (marked position on the package).

Break off the extraction buffer tube(b).


Drop all the extraction buffer into the specimen collection tube(a).

Take the Specimen Bag (c) and open it. Cough several times and clear your throat deeply. Spit sputum (about 2 ml) into the Specimen Bag (c). Remove the dropper and use it to collect samples.


Place 5 drops (approx. 200ul) of the saliva into the sample collection tube. Cover the sample collection tube, and put the rack.

TEST PROCEDURE

Place the test device on a clean and horizontal surface. Reverse the specimen collection tube, extrude 3 drops of the prepared specimen into the specimen well (S) of the test cassette and start the timer.

**Clinical Validation report of
SARS-CoV-2 Antigen Rapid Test
(COVID-19 Ag)**

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Jianfan HR

Final report date: 2021.02.02

Manufacturer: Hangzhou Immuno Biotech Co., Ltd.

I. Purposes

The SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) produced by Hangzhou Immuno Biotech Co., Ltd. was used to verify the feasibility of clinical evaluation and the reliability of test results compared with competitors.

II. Test design

1. Test plan selection and reasons

In vitro diagnostic reagents for testing and reference reagents were used to conduct comparative research tests on clinically suspected COVID-19 venous in human throat and nasal secretions, and it was proved that the in vitro diagnostic reagents used in the test can achieve the expected assistance in infection of the SARS-CoV-2.

2. Sample volume required

The total number of clinical trials of this product is not less than 200 cases. The samples is classified into the positive group and the negative group as per the test results of the reference product. Meanwhile, the samples shall be tested via the qualitative test strip tested and the reference one and then the test results of the product tested and the reference product shall be compared with statistical analysis being made.

3. Sample collection, processing and storage

Sample collection: Suitable for human throat, nasal secretions by the sterile swab, and saliva collected to the specimen bag in the morning.

Sample processing: Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing.

Sample storage: The assay should be performed immediately in 2 hours after the specimen preparation. If the assay could not be carried immediately, the prepared specimen should be kept no more than 24 hours at 2-8°C or 7 days at -20°C. Specimens should not be frozen and thawed repeatedly for more than two times. If specimens are to be shipped, they should be packed in compliance with federal regulations covering the transportation of etiologic agents.

4. In vitro diagnostic reagents and reference products for testing

4.1 Test in vitro diagnostic reagents

Name: SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) **Specification:** 20 tests/kit

LOT: 202008T01, 202009T01, 202012G15A

Storage Conditions: Store in a dry place at 2-30°C, protected from light. After opening the inner package, the test card will become invalid due to moisture absorption. Please use it within 1 hour.

Source: Hangzhou Immuno Biotech Co., Ltd.

4.2 Reference products

Name: Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit

Manufacturer: Shanghai ZJ Bio-Tech Co., Ltd., a commercial SARS-COV-2 kit approved by CFDA, is used as the "gold standard" reagent

III. method

1. Collect 200 cases of nasal swab, 100 cases of throat swab, and 100 saliva specimen from healthy people, and 107 cases of nasal swab, 110 cases of throat swab, and 123 cases of saliva specimen confirmed by the Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit.
2. The test group uses the in vitro diagnostic reagent products for testing, and the control group uses the "reference product" for testing.
3. The operation steps of the in vitro diagnostic reagents for the test are as follows. For details, please refer to the product instruction manual:

Step 1:

- A. For throat swab, insert the sterile swab into the throat. Gently scrape the secretions around the wall of pharynx for several times. Collect 0.5ml of assay buffer and place into a specimen collection tube. Insert the swab into the tube and squeeze the flexible tube to extrude the specimen from the head of the swab. Make the specimen resolved in the assay buffer sufficiently. Add the crystal tip onto the specimen collection tube. If the sample is stored refrigerated or frozen, remove the test sample and required reagents from the storage conditions and equilibrate to room temperature (15-30°C). After thawing, mix the samples thoroughly before testing.
- B. For nasal swab, insert the sterile swab into the deep nasal cavity. Gently rotate swab against wall of turbinate for several times. Make the swab wet as much as possible. Collect 0.5ml of assay buffer and place into a specimen collection tube. Insert the swab into the tube and squeeze the flexible tube to extrude the specimen from the head of the swab. Make the specimen resolved in the assay buffer sufficiently. Add the crystal tip onto the specimen collection tube.
- C. For saliva, Take a specimen collection container. Make a "Kruuuu" noise from the throat, to get out the saliva or sputum from the deep throat. Then spit saliva (about 1-2ml) into the container. Morning saliva is optimal for saliva collection. Do not brush the teeth, eat food or drink before collection the saliva specimen. Collect 0.5ml of assay buffer and place into a specimen collection tube. Suck the saliva from the container and place 5 drops (approx.200ul) of the saliva into the sample collection tube. Make the specimen resolved in the assay buffer sufficiently.

Step 2: When preparing for testing, open the aluminium foil bag from the tear. Remove the test card and lay it flat on a horizontal table.

Step 3: Label the sample number on the test card.

Step 4: Reverse the specimen collection tube, extrude 3 drops of the prepared specimen into the specimen well (S) of the test cassette and start the timer.

Step 5: Time counting and interpret the results within 10 minutes. Do not interpret the result after 15 minutes.

Note: The detection steps need to be completed under protection against infection.

IV. Statistical methods of statistical analysis of clinical research data

A Methods evaluating clinical performance

Whether various indexes can reach the standards of clinical evaluation shall be judged by calculating the consistency percentage of negative/positive and the total consistency percentage in the test results of the product tested and the reference product, to validate the accuracy and applicability of the product in clinical applications. The product tested shall be subject to tests through the sample of different types, with statistics on the results. Meanwhile, different types of sample of the subjects shall be subject to determination by the product tested synchronously, and then the determination results of both shall be compared. The test results recorded shall be subject to statistical analysis upon completion of determination of all clinical samples, to calculate the consistency percentage of negative/positive and the total consistency percentage. Afterwards, equivalence of both shall be evaluated as per these statistical indexes

B Statistical method

The products launched on the market shall be subject to comparative study and evaluation. Kappa inspection: each sample shall be tested with the product tested and the reference product respectively, and then the consistency in statistical results of these two inspection methods shall be compared through Kappa inspection.

The data shall be subject to Kappa inspection and analysis and the Kappa coefficient shall be calculated. Favorable consistency can be proven if Kappa is >0.8 . The consistency in test results of the product tested and the reference product is evaluated as per the evaluation standards.

V. Standards of clinical evaluation

The coincidence rate shall be calculated by comparing with the reference product whose marketing is approved. The product performance shall meet the following requirements.

- 1) Coincidence rate of negative: the sample whose test results are negative for both the product tested and the reference product and the proportion in the sample whose test results are negative for the reference product shall be more than 90%.
- 2) Coincidence rate of positive: the sample whose test results are positive for both the product tested and the reference product and the proportion in the sample whose test results are positive for the reference product shall be more than 90%.
- 3) Total coincidence rate: the sample whose test results are the same for the product tested and the reference product and its proportion in the total number of sample shall be more than 90%.

Method		Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit		Total Results
SARS-CoV-2 Antigen Rapid	Result	positive	negative	

Test (COVID-19 Ag)	positive	A	B	A+B
	negative	C	D	C+D
Total Results		A+C	B+D	A+B+C+D

Clinical sensitivity = $A/(A+C)*100\%$

Clinical specificity = $D/(B+D)*100\%$

Accuracy: $(A+D)/(A+B+C+D)*100\%$

If the coincidence rate of positive/negative can meet clinical requirements, two methods or Products are considered as equivalent; If the coincidence rate of positive/negative is greatly different, the clinical scheme should be re-designed.

VI. Provisions for amendments to clinical validation

In general, the clinical validation should not be changed. Any modification to the project during the test should be explained, and the time, reason, process of change, and whether there is a record of the change are explained in detail and its impact on the evaluation of the entire research result is explained.

VII. Results and Analysis of Clinical Tests

In total, 740 test samples are included for the unit and all test samples included are tested. Statistics on test results and those of the product tested are as follows:

Method		Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit		
		$C_T \leq 25$	$25 < C_T < 38$	$C_T \geq 38$
SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag)	positive	260	60	1
	negative	1	19	399
Coincidence Rate %		99.62%	75.95%	99.75%

*Consider PCR cycle threshold (C_T) value <38 as clinical positive.

Calculate the clinical sensitivity and specificity as below.

Method		Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit		Total Results
		positive	negative	
SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag)	Result	positive	negative	
	positive	320	1	321
	negative	20	399	419
Total Results		340	400	740

Clinical sensitivity = $320/340*100\% = 94.12\%$ (95%CI: 91.04%~96.21%)

Clinical specificity = $399/400*100\% = 99.75\%$ (95%CI: 98.45%~>99.99%)

Accuracy: $(320+399)/(320+20+1+399)*100\% = 97.16\%$ (95%CI: 95.68%~98.42%)

VIII. Discussion and Conclusions

1. Discussion

A Results of comparative analysis of the product tested and the reference product:

Test results of the nasal swab, throat swab and saliva sample of the product tested and the reference product: both the coincidence rate of negative/positive and the total coincidence rate are larger than 90%, indicating favourable consistency with the reference product.

B Statistical analysis results of the product tested for different types of clinical sample

While testing the SARS-CoV-2 antigen through the product tested for different types of clinical sample, the consistency percentages of negative/positive are 100.0% and the total consistency percentage is 100.0%, indicating favourable and complete consistency of two methods and equivalence of two such systems.

2. Test conclusions

By analysing the test results of the product tested and the reference product, the consistency percentage of negative/positive and the total consistency percentage are proven high. Moreover, according to the results of statistical analysis, there is no remarkable difference in test results of both, indicating favourable consistency in diagnosis and equivalence of two such systems. Meanwhile, the test results of the product tested for the serum and plasma sample of the same patient are completely identical. Therefore, such product is applicable to qualitative clinical analysis on the SARS-CoV-2 antigen in the nasal and throat secretions sample of humans, and can be used for auxiliary diagnosis of those suffering from pneumonia triggered by COVID-19.

IX. Quality control methods

On-site quality control

1) During the course of this study, clinical implementers appointed clinical inspectors to conduct regular on-site supervision visits to the research hospital. Through monitoring visits, it was found that all the contents of the research plan were strictly observed, and the correctness of the research data was also guaranteed. Participating researchers have undergone unified training, unified recording methods and judgment standards. The entire clinical trial process is conducted under strict operation, and the test content is complete and authentic. All observations and findings in the clinical trials have been verified and the data are reliable. The conclusions in the clinical trials are derived from the original data.

2) Quality control of clinical experiment process

During the evaluation, quality control was performed daily to ensure that the product was under control. Strict quality control is performed for each trial to ensure the quality of clinical trials.

X. Prediction of adverse events

Because the SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is an in vitro diagnostic reagent product, no direct contact with patients is required in clinical trials, no test report is provided to patients, and the test results are only used for comparative studies. It involves personal privacy, does not serve as a basis for auxiliary diagnosis, does not bring any risk to the subject, and does not cause adverse events.