PocRoc® SARS-CoV-2 ANTIGEN RAPID TEST KIT (Colloidal Gold)


PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) is intended for qualitative detection of SARS-CoV-2 nucleocapsid antigen from Swab specimen.

Specification:

- Target: nucleocapsid antigen
- Specimens: Nasal Swab (NS) / Oropharyngeal Swab(OP)
- Detection Time: 15~20mins
- Shelf Time: 24 months


Product Name	Specification	Storage temperature
PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)	1 test/box 5 tests/box 25 tests/box 50 tests/box	4°C-30°C

Lumigenex®


Lumigenex (Suzhou) Co., Ltd.

Sample preparation can take according to the operation steps.


- 1 The swab extractor bottle is pressed vertically downward to the extraction tube.
- 2 Put the swab specimen into the extraction tube, rotate the swab for about 10 seconds, and press the swab head against the tube wall to release the antigen in the swab.
- 3 Squeeze the swab over the head to remove the swab so as to remove as much liquid as possible from the swab.
- 4 Install the beater on the extraction tube.

Test procedure

Interpretation of results

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

**Lumigenex (Suzhou)
Co., Ltd.
Building C24
218 Xinghu Street
Suzhou Industrial Park
215123 Jiangsu
China**

has established and applies a quality management system for medical devices
for the following scope:

(see attachment for scope)

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2019-09-30
Certificate Registration No.: SX 60139875 0001
An audit was performed. Report No.: 15084639 004
This Certificate is valid until: 2021-12-19

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**Attachment to
Certificate**

Registration No.: SX 60139875 0001
Report No.: 15084639 004

Organization: Lumigenex (Suzhou)
Co., Ltd.
Building C24
218 Xinghu Street
Suzhou Industrial Park
215123 Jiangsu
China

Scope:

Design/development, Manufacture and Distribution of
In-vitro Diagnostic Medical Devices, including Time-Resolved
Fluorescence Immunoanalyzers and Test Kits; Dry Chemistry
Analyzers and Test Kits

Certification Body

Date: 2019-09-30

CERTIFICATE OF NOTIFICATION

This is to certify that, according to the European Council Directive 98/79/EC, Riomavix S.L. performed all notification duties and responsibilities as the European Authorized Representative:

MANUFACTURER: Lumigenex (Suzhou) Co., Ltd.

ADDRESS: Building C24, 218 Xing Hu Street, SIP, Suzhou, P.R. China 215123

The manufacturer has provided Riomavix S.L. with all the appropriate declaration according to the European Council Directive 98/79/EC including the Declaration of Conformity confirming that its in vitro diagnostic medical device, as stipulated here below, is fulfilling the essential requirements of the European Council Directive 98/79/EC.

IVD Devices:

PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Classification: Others

Where the manufacturer affix the CE mark to the device listed they must ensure that all the essential requirements of European Council Directive 98/79/EC are met.

The notification of abovementioned device has been completed by the European Authorized Representative in Spain. The Spain Competent Authority is notified of the manufacture's device and has allocated registration. The registration number is RPS/2506/2020

Executive Director

Issue date: 1/Nov/2020
Cert. No.: R20201101-1

Lumigenex

GJ-CE-2020DoC(Ver.02)

Declaration of Conformity

Manufacturer:

Lumigenex (Suzhou) CO., Ltd. Located at building C24, 218 Xing Hu Street, SIP, Suzhou, P.R. China 215123
TEL +86 (512) 80988088 FAX +86 (512) 80988096

European Representative:

Riomavix S.L.
Calle de Almansa 55, 1D, Madrid 28039 Spain

Products Name:

PocRoc[®] SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Product code: P23001, P23005, P23025, P23050.

Classification Under IVDD:

IVDD Other

Conformity assessment route: Annex III

We hereby declare under the sole responsibility of the manufacturer that the product mentioned above meets the provisions of the European Directive 98/79/EC for in vitro Diagnostic Medical Devices. All supporting documentations are retained at the premises of the manufacturer.

General applicable directive:

Directive 98/79/EC of European Parliament and of the Council of 27 October 1998 on in vitro Diagnostic Medical Devices.

Standards applied:

EN ISO 18113-1:2011, EN ISO 14971:2012, EN ISO 13485:2016, EN ISO 18113-2:2011, EN ISO 23640:2015, EN 13641:2002, EN ISO 15223-1:2016, EN ISO 17511:2000, EN 13612: 2002

First start of CE-MARC: Nov.1, 2020

Signature:

Mr. Eric Liu (Vice President)

A handwritten signature in purple ink that reads "Eric Liu".

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All data according to the transmission of the manufacturer, are only binding the information in the respective usage information.

Further information on the list provided by the BiArM as well as on the list and, if necessary, also removal from the list of underlying criteria can be found on our [website for antigen testing on SARS-CoV-2](#).

The following table shows the original tests with their trade name assigned by the manufacturer or European representative. An overview of the respective German distributors and their if necessary, different name can be found under the link in the column "German distributors".

The indication "Evaluation PEI" provides the corresponding overview of the comparative evaluation of the sensitivity of SARS-CoV-2 rapid antigen tests published on the website of the Paul Ehrlich Institute (PEI)(see [pei website](#)).

- "Yes" means that the test has already been evaluated with a positive result by the PEI.
- "No" means that no corresponding test results are available yet.

In the event of a negative evaluation by the PEI, the BiArM deletes the corresponding test with all associated distributors from its list.

		Manufacturer			European Representative			Sensitivity		Specificity			
Test ID	Manufacturer's trade name / Europ. Agent	Name	City	Co...	Name	City	Co...	German distribut...	Test locat...	%	95% confiden... Interval	%	95% confidence Interval
AT45520	PocRoc® SARS-CoV-2 Antigen Quick Detection Kit (Gold Colloidal)	Lumigenex (Suzhou) Co., Ltd.	Suzhou	Cn	Riomavix S.L.	Madrid	It	De...	POC (without device)	93,33	86,27 - 97,05	99,16	96,68 - 99,85

Actions

Los

Lumigenex

Reset

Search for 'Lumigenex'

COVID-19 In Vitro Diagnostic Devices and Test Methods Database

Home > COVID-19 In Vitro Diagnostic Medical Devices > COVID-19 In Vitro Diagnostic Medical Device - detail

COVID-19 In Vitro Diagnostic Medical Device - detail

PocRoc®SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

Manufactured by Lumigenex (Suzhou) Co., Ltd, China - <http://www.lumigenex.com/>

Device identification number	2128
CE Marking	✔ Yes
HSC common list	✔ Yes
HSC mutual recognition	✔ Yes
Format	Manual, Near POC / POC
Physical Support	Cassette, Lateral flow
Target	Antigen