Exdia COVID-19 Ag

One Step Immunoassay for Exdia COVID-19 Ag For *in vitro* Diagnostic Use

One-step qualitative time-resolved fluorescence immuno-chromatography assay for the detection of COVID-19 antigen in human nasopharyngeal secretion

Manufactured by Precision Biosensor Inc.

1. INTENDED USE

Exdia COVID-19 Ag is a time-resolved fluorescence immunoassay that is used with the Exdia TRF analyzer intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider. The Exdia COVID-19 Ag test identifies SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral antigens but clinical correlation with patient history and other diagnostic information are necessary to determine the status of infection. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Negative results should be treated as presumptive, do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management. The Exdia COVID-19 Ag is intended for use by trained clinical laboratory personnel specifically instructed and trained *in vitro* diagnostic procedures and individuals similarly trained in point of care settings.

2. SUMMARY AND EXPLANATION

SARS-CoV-2, which causes the disease COVID-19, was first identified in Wuhan, Hubei Province, China, December 2019. The WHO declared that COVID-19 was pandemic on March 11, 2020 and human infection has spread globally with hundreds of thousands of confirmed infections and deaths. The median incubation time is estimated to be approximately 5 days with symptoms estimated to be present 12 days after the infection. The symptoms of COVID-19 are similar to other viral respiratory diseases including fever, cough, shortness of breath.

3. PRINCIPLE OF THE TEST

The Exdia COVID-19 Ag employs the time-resolved fluorescence technology in a sandwich assay design that is used with Exdia TRF analyzer to detect nucleocapsid protein antigen from SARS-CoV-2. This allows for the detection of SARS-CoV and SARS-CoV-2. The test detects, but does not differentiate, between the two viruses. The patient sample is placed in the Reagent tube, where the virus particles in the sample are disrupted, exposing internal viral nucleoproteins. After disruption, the sample is dispensed into the sample well of the test cassette. The present SARS antigens form the immunocomplexes with both the biotinylated and fluorescence nanoparticle-conjugated anti-SARS antibodies during the migration on the cassette membrane. These immunocomplexes bind to streptavidin-immobilized test line resulting in the fluorescent signal. Exdia TRF analyzer will measure the fluorescent signal and display the test results (Positive, Negative, or Invalid).

4. REAGENT AND MATERIALS

Provided

- 20 Test cassettes sealed in pouch with desiccant
- 20 Sterilized nasopharyngeal swabs for sample collection
- 20 Reagent tubes with extraction buffer (0.3 mL/tube)
- 20 Filter caps
- 2 Control swabs
 - 1 COVID-19 Positive control swab
 - 1 COVID-19 Negative control swab
- 1 QR card for calibration (lot-specific)
- 1 Instructions for Use

Required but not provided

- Timer
- Exdia TRF analyzer
- Exdia IQC cassette

5. STORAGE AND STABILITY

The test cassette should be stored at 2°C ~ 30°C in the original sealed pouch for the duration of shelf life.

6. PRECAUTIONS

- For *in vitro* diagnostic and professional use only.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be established.
- To avoid cross contamination, use a fresh nasopharyngeal swab for each clinical sample tested.
- Do not use test kit if the pouch is damaged or improperly sealed.
- Do not use test kit beyond expiration date.
- Use only the swabs provided. Other swabs may not work properly.
- Extraction buffer contains slightly caustic chemicals. Avoid contact with eyes, sensitive mucous membranes. If the reagent comes in contact with skin or eyes, flush with a large volume of water.
- To avoid contact with contagious agents, wear disposable gloves while handling kit reagents or specimens and thoroughly wash hands afterwards.


7. SPECIMEN COLLECTION AND PREPARATION

- To collect the nasopharyngeal swab sample, carefully insert the swab into the nostril that presents the most secretion under visual inspection. Keep the swab near the septum floor of the nose while gently pushing the swab into the posterior nasopharynx. Rotate the swab several times then remove it from the nasopharynx.
- The samples should be collected under the standard laboratory conditions.
- Samples should be tested immediately after collection.
- Refrigerated or frozen specimen should be brought to room temperature (20°C ~ 30°C) and be homogeneous prior to testing.

※ Use of VTM (viral transport media) or UTM (universal transport media) may result in decreased test sensitivity due to dilution, and the direct testing of swab from the patient is highly recommended.

8. TEST PROCEDURE AND PROTOCOL

- Collect specimen according to the instructions in “SPECIMEN COLLECTION AND PREPARATION”.
- Test cassette and sample should be brought to room temperature (20°C ~ 30°C) prior to testing.
- Remove the test cassette from the sealed pouch immediately before use. Label the cassette with patient or control identification.
- Open the Reagent tube and discard the tube cap.
- Insert the collected swab into the Reagent tube and rotate the swab against the inner wall of the tube for several times (more than five times). Plunge the swab up and down in the extraction buffer while squeezing the sides of the tube for 15 seconds.
- Remove the swab while squeezing the sides of the tube to extract the maximum amount of liquid from the swab. Properly discard the swab.
- Close the Reagent tube containing the sample firmly using the Filter cap.
- Invert the Reagent tube and hold the tube vertically (approximately 1 inch above sample well of the test cassette). While holding middle part of the tube, squeeze gently to dispense two drops of the sample into the sample well of the test cassette.
- The test cassette should be analyzed by the Exdia TRF analyzer following the instructions in the manual. The test cassette will be analyzed at the different preset times after applying the samples. If the test result is strong positive, the analyzer will immediately determine the result as positive. If the test result is positive at near cutoff level or negative, the result will be determined at 20 min after applying the sample. In skip mode, read the results at 20 minutes (19.5 min ~ 20.5 min).

< Regarding VTM or UTM >

If the previously collected samples in VTM or UTM must be used for Exdia COVID-19 Ag, minimal dilution of the sample is recommended. Mix previously collected samples in VTM/UTM to the extraction buffer at a one to one ratio (mixing volume is 300 µL each). After mixing for 1 minute, close the Filter cap, and follow the subsequent test procedure.

9. INTERPRETATION OF RESULTS

Exdia COVID-19 Ag results are interpreted as qualitative method by Exdia TRF analyzer. The qualitative results are displayed as Negative or Positive or Invalid based on the measured fluorescence intensity of test and control line. If the test result is Invalid, a new test should be performed with a new test cassette.

10. LIMITATIONS

- The test is for professional and *in vitro* diagnostic use only.
- Use of VTM or UTM may result in decreased test sensitivity due to dilution, and direct testing of swab is highly recommended.
- As with all antigen tests, performance may decrease with the number of days after symptom onset.
- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the results of a single test. The test result should be used in conjunction with other clinical information such as clinical signs and symptoms and other test results to diagnose SARS-CoV-2 infection. Confirmation of test results should only be made by a physician along with clinical symptoms and laboratory findings.
- This test detects both viable (live) and non-viable SARS-CoV and SARS-CoV-2 virus. Test performance depends on the amount of virus (antigen) in the sample and may not correlate with viral culture results performed on the same sample.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly.
- Failure to follow the test procedure may adversely affect test performance and/or invalidate the test result.
- Test results must be evaluated in conjunction with other clinical data available to the physician.
- Positive test results do not rule out co-infection with other pathogens.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative results, from patients with symptom onset beyond five days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.

11. QUALITY CONTROL

The presence of fluorescence band in the control area of the window acts as an internal control to ensure an adequate volume of sample has been added. In the absence of this control band, the associated test result is invalid and the sample must be retested. Good laboratory practice recommends quality control to ensure proper test performance. External positive and negative control swabs are supplied in the test. The COVID-19 positive control swab contains SARS-CoV-2 nucleoprotein antigen. The positive and negative control swabs must be regularly run. If the quality control procedures in your laboratory require more frequent use of controls to verify the test results, follow your laboratory-specific procedures. The regular testing of Exdia TRF analyzer using Exdia IQC cassette is recommended depending on the laboratory practice. If the test result is questionable for any reason, contact the customer support.

12. PERFORMANCE CHARACTERISTICS

12-1. Clinical performance study

The clinical performance of the Exdia COVID-19 Ag was established with a study using ninety-nine (99) previously characterized frozen nasopharyngeal swabs originally collected in transport media. The Exdia COVID-19 Ag demonstrated the positive percent agreement (PPA) of 93.9% and negative percent agreement (NPA) of 98.0% when compared to an EUA RT-PCR method (Allplex™ 2019-nCoV Assay, Seegene).

		RT-PCR (Allplex™ 2019-nCoV Assay)		Total
		POS	NEG	
Exdia COVID-19 Ag	POS	46	1	47
	NEG	3	49	52
	Total	49	50	99

- Positive Percent Agreement (PPA, %) = (46/49) x 100 = 93.9% (Wilson 95% CI: 83.5-97.9)
- Negative Percent Agreement (NPA, %) = (49/50) x 100 = 98.0% (Wilson 95% CI: 89.5-99.6)

Based on the positive percent agreement study, the results were further evaluated according to RT-PCR threshold cycle (Ct) values.

PCR Ct (RdRp gene)	Sample number	Number of positive of Exdia COVID-19 Ag (PPA, %)	Wilson 95% CI
Ct ≤ 20	6	6 (100)	61.0-100
20 < Ct ≤ 30	34	33 (97.1)	85.1-99.5
Ct > 30	9	6 (66.7)	35.4-87.9

12-2. Supportive contrived performance study

Due to the limited availability of direct swabs, spiking of heat-inactivated virus into negative nasopharyngeal swab was performed to simulate the direct swab testing procedure and performance. Thirty nasopharyngeal swabs were spiked with approximately 50 µL of a virus dilution. Each ten (10) specimens were respectively spiked with 1x LoD (7.55 x 10² TCID₅₀/mL), 2x LoD (1.5 x 10³ TCID₅₀/mL), and 3x LoD (2.27 x 10³ TCID₅₀/mL) of virus. The spiked swabs were added to the Exdia COVID-19 extraction solution concurrently to thirty nasopharyngeal swabs from thirty different individuals. Thirty (30) negative nasopharyngeal swabs from thirty different individuals were also tested. The table below presents the data for this study:

Sample type	Sample number	Number of Positives/Number of Tested
Negative (unspiked)	30	0/30
1x LoD	10	10/10
2x LoD	10	10/10
3x LoD	10	10/10

- Positive Percent Agreement (PPA, %) = (30/30) x 100 = 100% (Wilson 95% CI: 88.6-100)
- Negative Percent Agreement (NPA, %) = (30/30) x 100 = 100% (Wilson 95% CI: 88.6-100)

12-3. Limit of Detection (LoD)

The Limit of Detection (LOD) of Exdia COVID-19 Ag was determined using serial dilution of heat-inactivated SARS-CoV-2 (Zepmetrix, Isolate: USA-WA1/2020). The material was supplied frozen at a concentration of 1.51 x 10⁶ TCID₅₀/mL. The study to determine the Exdia COVID-19 Ag LOD was designed to reflect the assay procedure when using direct swab. Nasopharyngeal swab was spiked with approximately 50 µL of virus dilution in saline. The spiked nasopharyngeal swab was added to the Reagent tube with extraction buffer according to the test procedure.

The test consisted of two steps for LoD determination:

Tentative LoD screening

10-fold serial dilution of the heat-inactivated virus were made, and they were tested in 5 replicates. The concentration demonstrating 5 of 5 positives were chosen to be 7.55 x 10² TCID₅₀/mL. And 2-fold serial dilution was further performed and tested to estimate LoD.

LoD confirmation

The sample dilution with a final concentration of 7.55 x 10² TCID₅₀/mL was tested further. Twenty (20) of twenty (20) results were found to be positive. Based on this testing, the LoD was concluded to be 7.55 x 10² TCID₅₀/mL.

Dilution	SARS-CoV-2 (Heat-inactivated, USA-WA1/2020)								
	1/1	1/10	1/100	1/1,000	1/2,000	1/4,000	1/8,000	1/100,000	1/151,000
Concentration (TCID ₅₀ /mL)	1.51 x 10 ⁶	1.51 x 10 ⁵	1.51 x 10 ⁴	1.51 x 10 ³	7.55 x 10 ²	3.78 x 10 ²	1.89 x 10 ²	1.51 x 10 ¹	1.00 x 10 ¹
Positive call rate	100% (5/5)	100% (5/5)	100% (5/5)	100% (5/5)	100% (20/20)	5% (1/20)	0% (0/20)	0% (0/20)	0% (0/20)
LoD	7.55 x 10 ² TCID ₅₀ /mL								

12-4. Cross-reactivity

Cross-reactivity and potential interference of the Exdia COVID-19 Ag was evaluated by testing various microorganism or virus samples. Each microorganism and virus were tested in triplicate in the absence or presence of the heat-inactivated SARS-CoV-2 virus 1.51 x 10³ TCID₅₀/mL.

No.	Virus/Bacteria	Strain	Sample type	Concentration	Cross-reactivity results*	Interference results*
1	Adenovirus A	Type 18 (KUMC-4)	Isolate	1.2 x 10 ⁶ pfu/mL	NEG	POS
2	Adenovirus B	Type 11 (KUMC-63)	Isolate	3.0 x 10 ⁵ pfu/mL	NEG	POS
3	Adenovirus C	Type 5 (KUMC-61)	Isolate	4.0 x 10 ⁶ pfu/mL	NEG	POS
4	Adenovirus D	Type 23 (KUMC-5)	Isolate	1.2 x 10 ⁶ pfu/mL	NEG	POS
5	Adenovirus E	Type 4 (KUMC-60)	Isolate	2.0 x 10 ⁵ pfu/mL	NEG	POS
6	Adenovirus F	Type 40 (KUMC-6)	Isolate	3.0 x 10 ⁵ pfu/mL	NEG	POS
7	Alpha Coronavirus	229E (KUMC-9)	Isolate	8.0 x 10 ⁵ pfu/mL	NEG	POS
8	Alpha Coronavirus	NL63	Isolate	1.0 x 10 ⁵ pfu/mL	NEG	POS
9	Beta Coronavirus	OC43 (KUMC-8)	Isolate	9.8 x 10 ⁵ pfu/mL	NEG	POS

10	MERS Coronavirus	EMC/2012 (BEI NR-50549)	Isolate	1.27 x 10 ⁵ TCID ₅₀ /mL	NEG	POS
11	<i>Mycoplasma pneumoniae</i>	M129-B7	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
12	<i>Streptococcus pyogenes Rosenbach</i>	Bruno (CIP 104226)	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
13	Influenza A H3N2	KUMC-32	Isolate	4.0 x 10 ⁵ pfu/mL	NEG	POS
14	Influenza A H1N1	KUMC-33	Isolate	1.2 x 10 ⁵ pfu/mL	NEG	POS
15	Influenza B	KUMC-34	Isolate	4.9 x 10 ⁵ pfu/mL	NEG	POS
16	Human Parainfluenza virus 1	KUMC-64	Isolate	2.0 x 10 ⁵ pfu/mL	NEG	POS
17	Human Parainfluenza virus 2	KUMC-65	Isolate	2.0 x 10 ⁵ pfu/mL	NEG	POS
18	Human Parainfluenza virus 3	KUMC-67	Isolate	4.0 x 10 ⁶ pfu/mL	NEG	POS
19	Human Parainfluenza virus 4a	KUMC-69	Isolate	4.6 x 10 ⁷ pfu/mL	NEG	POS
20	Human Enterovirus A	Type 71 (KUMC-56)	Isolate	0.7 x 10 ⁵ pfu/mL	NEG	POS
21	Human Enterovirus B	Coxsackievirus B3 (KUMC-15)	Isolate	8.6 x 10 ⁶ pfu/mL	NEG	POS
22	Human Enterovirus C	Poliovirus 1 (Sabin)	Isolate	6.0 x 10 ⁶ pfu/mL	NEG	POS
23	Human Enterovirus D	Type 70 (KUMC-55)	Isolate	8.8 x 10 ⁶ pfu/mL	NEG	POS
24	Human Metapneumovirus	KUMC-87	Isolate	1.4 x 10 ⁵ pfu/mL	NEG	POS
25	Human Respiratory Syncytial Virus	KUMC-41	Isolate	8.0 x 10 ⁵ pfu/mL	NEG	POS
26	Human Rhinovirus A	Type 1B (KUMC-81)	Isolate	1.2 x 10 ⁵ pfu/mL	NEG	POS
27	Human Rhinovirus B	Type 42 (KUMC-80)	Isolate	4.2 x 10 ⁵ pfu/mL	NEG	POS
28	<i>Candida albicans</i>	4454M	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
29	<i>Chlamydomonas pneumoniae</i>	CM-1	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
30	<i>Legionella pneumophila</i>	Bloomington-2	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
31	<i>Streptococcus pneumoniae</i>	262 (CIP 104340)	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS
32	<i>Bordetella pertussis</i>	MN2531	Isolate	5.0 x 10 ⁴ cells/mL	NEG	POS

* Testing was performed in triplicate.

12-5. Interference substance study

Potentially interfering substances were spiked into SARS-CoV-2 positive or negative samples. The substances at the following level do not interfere with the performance of the Exdia COVID-19 Ag. No false positive or negative result was found with following:

Solvent	Substance	Final conc.	Cross-reactivity results*	Interference results*
DMSO	Acetaminophen	10 mg/mL	NEG	POS
	Bilirubin	15 mg/mL	NEG	POS
PBS pH 7.2	Tobramycin	51.4 µM	NEG	POS
Methanol	Dexamethasone	1.53 µM	NEG	POS
DW	Oseltamivir	5 mg/mL	NEG	POS
	Human serum albumin	50 mg/mL	NEG	POS
	Glucose	1.2 mg/mL	NEG	POS
NaOH	Mucin	1 mg/mL	NEG	POS
Ethanol	Benzocaine	1 mg/mL	NEG	POS
PBS pH 7.2	Whole blood	1%	NEG	POS

* Testing was performed in triplicate.

12.6 Hook effect

The highest concentration of heat-inactivated SARS-CoV-2 stock available (1.51 x 10⁶ TCID₅₀/mL) was tested. There was no Hook effect observed.

13. REFERENCES

1. Baker, S., Frias, L., and Bendix, A. Coronavirus live updates: More than 92,000 people have been infected and at least 3,100 have died. The US has reported 6 deaths. Here's everything we know. Business Insider. March 03, 2020.
2. <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/events-as-they-happen>

3. Lauer, S.A., et al. The incubation period of Coronavirus disease 2019 (COVID-19) from publicly reported confirmed cases: estimation and application, Ann Intern Med. 2020

For more information or any questions about this product, please contact customer support at: support@precision-bio.com

 **Precision Biosensor Inc.**
306, Techno 2-ro, Yuseong-gu,
Daejeon, 34036
Republic of Korea
Tel: +82-42-867-6300
Fax: +82-42-867-6302
www.precision-bio.com

 **mdi Europa GmbH**
Langenhagener Str. 71, 30855
Langenhagen, Germany
Tel: +49 511 39 08 95 30
e-mail: info@mdi-europa.com
www.mdi-europa.com

Document Number: PB-FGRUM-01-EN
Revision Number: 01
Date of effective: Oct 27, 2020

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30	<i>Legionella pneumophila</i>	Bloomington-2	Isolat	$5,0 \times 10^4$ cells/mL	NEG	POS
31	<i>Streptococcus pneumoniae</i>	262 (CIP 104340)	Isolat	$5,0 \times 10^4$ cells/mL	NEG	POS
32	<i>Bordetella pertussis</i>	MN2531	Isolat	$5,0 \times 10^4$ cells/mL	NEG	POS

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13. Literatur

- Baker, S., Frias, L., and Bendix, A. Coronavirus live updates: More than 92,000 people have been infected and at least 3,100 have died. The US has reported 6 deaths. Here's everything we know. Business Insider. March 03, 2020.
- <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/events-as-they-happen>
- Lauer, S.A., et. al. The incubation period of Coronavirus disease 2019 (COVID-19) from publicly reported confirmed cases: estimation and application, Ann Intern Med. 2020

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