# SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) Package Insert

REF:HRK-6820

**A RAPID TEST FOR THE QUALITATIVE DETECTION OF SARS-COV-2 ANTIGENS IN HUMAN SALIVA.**

## INTENDED USE

SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) is an *in vitro* diagnostic test for the qualitative detection Nucleoprotein and spike glycoprotein of Coronavirus Disease 2019 in human Oropharyngeal saliva, using the rapid immunochromatographic method as an aid in the diagnosis of SARS-CoV-2 infections. The identification is based on the monoclonal antibodies specific for the SARS-CoV-2 antigen. It will provide information for clinical doctors to prescribe correct medications.

## SUMMARY

The SARS-CoV-2 belong to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the SARS-CoV-2 are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Severe acute respiratory syndrome - coronavirus- 2 (SARS-CoV-2) is an enveloped, non-segmented Positive sense rna virus. It is the cause of the Coronavirus-0 disease (COVID-19) common to humans is contagious. SARS-CoV-2 has several structural proteins, including spike (S), envelope (E), membrane (M) and nucleocapsid(N).

At present, there are many variants of the SARS-CoV-2, and the N501Y mutation and its approximate variants have attracted attention because their mutation position is located in the spike glycoprotein's receptor-binding domain of the virus, thereby changing the virus infected efficiency. In silico analysis demonstrated that the N501Y mutation did not alter the primary and tertiary protein structure of the spike protein RBD domain. Therefore, its antigenicity remains unchanged.

## PRINCIPLE

SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to SARS-CoV-2

The test strip is composed of the following three parts, namely sample pad, reagent pad and reaction membrane. The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against SARS-CoV-2; the reaction membrane contains the secondary antibodies for SARS-CoV-2, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

When the saliva sample is received by the test, the conjugated solution from the reagent pad gets dissolved and migrates along with the saliva. When the SARS-CoV-2 is present in the saliva sample, a complex is formed between the anti-SARS-CoV-2 conjugate and the virus will be caught / detected by the specific anti- SARS-CoV-2 monoclonal coated on the T region. Whether the sample contains the virus or not, the solution continues to migrate to encounter another reagent (an anti-mouse IgG antibody) that binds the remaining conjugates, thereby producing a red line on the region C.

SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) can detect both the SARS-Cov-2 nucleoprotein as well as the SARS-Cov-2 spike protein. By ELISA, we determined that the antibody we use binds to amino acids 511-531 of the SARS Cov-2 spike protein.

The detectability of genetic SARS-CoV-2 variants was tested by examining the sensitivity toward recombinant SARS-Cov-2 spike proteins (319 to 541aa). In these tests, SARS-CoV-2 Antigen Saliva Rapid Test Kit achieved the same values when detecting the B.1.1.7 (UK) and B.1.351 (SA) variants as when detecting the standard variant.

## REAGENTS

The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against SARS-CoV-2; the reaction membrane contains the secondary antibodies for SARS-CoV-2, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

## PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- Do not use after the expiration date.
- Ensure foil pouch containing test device is not damaged before opening for use.
- Perform test at room temperature 15 to 30°C.
- Wear gloves when hanging the samples, avoid touching the reagent membrane and sample window.
- All samples and used accessories should be treated as infectious and discarded according to local regulations.

- Avoid using bloody samples.

## STORAGE AND STABILITY

SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) at room temperature or refrigerated (2-30°C). Do not freeze. All reagents are stable until the expiration dates marked on their outer packaging and buffer vial.

## SPECIMEN COLLECTION AND PREPARATION

### 1. Specimen collection:

The oral fluid specimen should be collected using the collection tools provided with the kit. Follow the detailed Directions for Use below. No other collection tools should be used with this assay. Oral fluid collected at any time of the day may be used.

### 2. Specimen preparation:

When the saliva is collected, follow the direction to prepare the specimen with buffer provided with the kit.

## MATERIALS

### Materials provided

- 20\*Test Cassette
- 20\*Saliva Collector
- 20\*Extraction buffer
- 1\*Package Insert
- 20\*Extraction Tube
- Tube stand\*

\*The 20-test package contains the tube stand, the 1-test and 5-test package use the test box itself as tube stand.



### Materials required but not provided

- Timer

## DIRECTIONS FOR USE

**Allow the test device, specimen, extraction buffer to equilibrate to room temperature (15-30°C) prior to testing. Do not place anything in the mouth including food, drink, gum, tobacco, water and mouthwash products for at least 10 minutes prior to collection of oral fluid specimen.**

1. Screw the saliva collector onto the extraction tube.
2. Place the tip of the tongue against the mucous membrane of the oral cavity to enrich the saliva, and then spit the saliva into a saliva collector, reach the 300ul line of extraction tube (Please don't exceed the line)
3. Hold the saliva collection tube vertically and screw the saliva collector to remove it from the tube. Dispose of the funnel.
4. Take out a bottle of extraction buffer, remove the extraction buffer bottle cap, add all the extraction buffer into the extraction tube.
5. Close the tube tightly with the filter cap. Mix the contents by turning the tube ten times.
6. Remove the test device from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch. Put the test device on a clean and flat surface.
7. Transfer 3 drops of sample into the sample well of test device vertically, start the timer.
8. Read the result at 10-20 minutes. Don't interpret the result after 20 minutes.

## INTERPRETATION OF RESULTS

(Please refer to the illustration above)

**POSITIVE:** Two red lines appear. One red line appears in the control region(C), and one red line

in the test region(T). The shade of color may vary, but it should be considered positive whenever there is even a faint line.

**NEGATIVE:** Only one red line appears in the control region(C), and no line in the test region(T). The negative result indicates that there are no SARS-CoV-2 particles in the sample or the number of viral particles is below the detectable range.

**INVALID:** No red line appears in the control region(C). The test is invalid even if there is a line on test region(T). Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the test procedure and repeat the test using a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

## LIMITATIONS

• SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) is an acute-phase screening test for qualitative detection. Sample collected may contain antigen concentration below the reagent's sensitivity threshold, so a negative test result does not exclude infection with SARS-CoV-2.

• SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) detects viable and non-viable SARS-CoV-2 antigen. Test performance depends on antigen load in the sample and may not correlate with cell culture performed on the same sample. A positive test does not rule out the possibility that other pathogens may be present, therefore, the results must be compared with all other available clinical and laboratory information to make an accurate diagnosis.

• A negative test result may occur if the level of extracted antigen in a specimen is below the sensitivity of the test or if poor quality specimen is obtained.

• Performance of the test has not been established for monitoring antiviral treatment of SARS-CoV-2.

• Positive test results do not rule out co-infections with other pathogens.

• Negative test results are not intended to rule in other coronavirus infection except the SARS-Cov-2.

• Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children List.

• The concentration of virus in saliva is greatly affected by factors such as meals, diet, smoking, breath fresheners, etc. Therefore, please strictly follow this manual before collecting samples. A negative result may occur if the concentration of antigen in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-Cov-2 infection, and should be confirmed by viral culture or PCR.

## PERFORMANCE CHARACTERISTICS

### Clinical Evaluation

Clinical evaluation was performed to compare the results obtained by SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) and PCR test result.

The performance has been determined using 290 Oropharyngeal saliva from individual symptomatic patients with Suspected COVID-19. Oropharyngeal saliva was made refer to the instructions for use of the specimen collection and preparation. All samples were selected and then tested sequentially in a blinded fashion. The performance of the test device was matched with the results of a compared to commercial RT-PCR reagent. The test device showed a sensitivity of 97.64% and a specificity of 100%. The results were summarized below:

Table: SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) vs. PCR

Method	SARS-CoV-2 Nucleic Acid Test Kit (RT-PCR)		Total Results	
	Results	Positive		Negative
SARS-CoV-2 Antigen Saliva Rapid Test Kit	Positive	124	0	124
	Negative	3	163	166
<b>Total Results</b>		127	163	290

Clinical sensitivity = 124/127= 97.64 % (95%CI\*:93.25% to 99.51%)

Clinical specificity =163/163=100% (95%CI\*:97.76% to >100.00%)

Accuracy: (124+163)/(124+3+163) \*100%=98.97% (95%CI\* 97.01% to 99.79%)

\*Confidence Interval

### Limit of Detection (LOD)

2019-nCoV Strain Tested	Huaree Tech product				
Stock 2019-nCoV Concentration	1 X 10 <sup>2</sup> TCID <sub>50</sub> /mL				
Dilution	1/100	1/200	1/400	1/800	1/1600
Concentration in Dilution tested (TCID <sub>50</sub> /ml)	1X10 <sup>3</sup>	5X10 <sup>2</sup>	2.5X10 <sup>2</sup>	1.25X10 <sup>2</sup>	62.5
Call rates of 20 replicates near cut-off	100(20/20)	100(20/20)	100(20/20)	95(19/20)	10(2/20)
Limit of detection (LOD) per Virus Strain	1.25 X 10 <sup>2</sup> TCID <sub>50</sub> /mL				

### Cross Reaction

The Cross reaction study conducted By using SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography) to test sample in which add each concentration of pathogens listed in the following table to the specimens prepared negative samples and 3XLOD positive samples. The results show that the pathogenic listed in the following table have no cross-reaction that would create false positive or negative results for SARS-Cov-2 antigen.

Virus/Bacteria/Parasite	Strain	Concentration
MERS-coronavirus	N/A	72µg/mL
	Type 1	1.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Adenovirus	Type 3	7.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 5	4.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 7	1.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL

	Type 8	1.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 11	2.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 18	2.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 23	6.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 55	1.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza A	H1N1 Denver	3.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	H1N1 WS/33	2.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	H1N1 A/Mal/302/54	1.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	H1N1 New Caledonia	7.6 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	H3N2 A/Hong Kong/8/68	4.6 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Influenza B	Nevada/03/2011	1.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	B/Lee/40	8.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
	B/Taiwan/2/62	4.0 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Respiratory syncytial virus	N/A	2.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Legionella pneumophila	Bloomington-2	1 x 10 <sup>5</sup> PFU/mL
	Los Angeles-1	1 x 10 <sup>5</sup> PFU/mL
	82A3105	1 x 10 <sup>5</sup> PFU/mL
Rhinovirus A16	N/A	1.5 x 10 <sup>6</sup> TCID <sub>50</sub> /mL
Mycobacterium tuberculosis	K	1 x 10 <sup>5</sup> PFU/mL
	Erdman	1 x 10 <sup>5</sup> PFU/mL
	HN878	1 x 10 <sup>5</sup> PFU/mL
	CDC1551	1 x 10 <sup>5</sup> PFU/mL
	H37Rv	1 x 10 <sup>5</sup> PFU/mL
Streptococcus pneumonia	4752-98 [Maryland (D1)6B-17]	1 x 10 <sup>5</sup> PFU/mL
	178 [Poland 23F-16]	1 x 10 <sup>5</sup> PFU/mL
	262 [CIP 104340]	1 x 10 <sup>5</sup> PFU/mL
Streptococcus pyrogens	Slovakia 14-10 [29055]	1 x 10 <sup>5</sup> PFU/mL
	Typing strain T1 [NCIB 11841, SF 130]	1 x 10 <sup>5</sup> PFU/ml
Mycoplasma pneumoniae	Mutant 22	1 x 10 <sup>5</sup> PFU/ml
	FHstrain of EatonAgent [NCTC 10119]	1 x 10 <sup>5</sup> PFU/ml
	36M129-B7	1 x 10 <sup>5</sup> PFU/ml
Coronavirus	229E	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	OC43	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	NL63	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	HKU1	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
Human etapneumovirus (hMPV) 3 Type B1	Peru2-2002	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
Human Metapneumovirus (hMPV) 16 Type A1	IA10-2003	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
Parainfluenza virus	Type 1	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 2	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 3	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL
	Type 4A	1.5 x10 <sup>6</sup> TCID <sub>50</sub> /mL

#### Interfering Substances Reaction



When tested using SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography), there was no interference between the device reagents and the Potential interference substances listed in below table that would create false positive or negative results for SARS-Cov-2 antigen.

Substance	Concentration	Substance	Concentration
Mucin	100µg/mL	Acetylsalicylic acid	3.0 mM
Whole Blood	5% (v/v)	Ibuprofen	2.5 mM
Biotin	100µg/mL	Mupirocin	10 mg/mL
Neo-Synephrine (Phenylephrine)	5%(v/v)	Tobramycin	10µg/mL
Afrin Nasal Spray (Oxymetazoline)	5%(v/v)	Erythromycin	50µM
Saline Nasal Spray	5%(v/v)	Ciprofloxacin	50µM
Homeopathic	5%(v/v)	Ceftriaxone	110mg/mL
Sodium Cromoglycate	10 mg/mL	Meropenem	3.7µg/mL
Olopatadine Hydrochloride	10 mg/mL	Tobramycin	100µg/mL
Zanamivir	5 mg/mL	Histamine Hydrochloride	100µg/mL
Oseltamivir	10 mg/mL	Peramivir	1mmol/mL
Artemether-lumefantrine	50µM	Flunisolide	100µg/mL
Doxycycline hyclate	50µM	Budesonide	0.64nmol/ L
Quinine	150µM	Fluticasone	0.3ng/mL
Lamivudine	1 mg/mL	Lopinavir	6µg/mL
Ribavirin	1 mg/mL	Ritonavir	8.2mg/mL
Daclatasvir	1 mg/mL	Abidor	417.8ng/mL
Acetaminophen	150µM	Pooled human mouth wash	N/A

#### Dose Hook Effect


Test SARS-CoV-2 wild strain culture medium (concentration 1X10<sup>7</sup> TCID<sub>50</sub>/ml) and multiple dilution sample with SARS-CoV-2 Antigen Saliva Rapid Test Kit (Immunochromatography), no hook effect was found.

#### SYMBOL

Symbol	Meaning	Symbol	Meaning
	In vitro diagnostic medical device		Storage temperature limit
	Manufacturer		Authorized representative in the European Community
	Date of Manufacture		Use by date
	Do not reuse		Consult instruction for use
	Batch code		Meet the requirements of EC Directive 98/79/EC

**SHENZHEN HUAREE TECHNOLOGY CO.,LTD**  
 3th Building, Huahang Technology Industrial Park,  
 No.19 Lanjin 4th road, Heping Community, Pingshan  
 District, Shenzhen 518118, China  
**Website: www.huareetech.com**  
**Email: info@huareetech**

**SUNGO Europe B.V.**  
 Olympisch Stadion 24, 1076DE  
 Amsterdam, Netherlands  
**Email: ec.rep@sungogroup.com**  
**Tel: +31(0)2021 11106**

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