Deepblue COVID-19 Antigentest zur Laienanwendung

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Hersteller

ANHUI DEEPLBLUE MEDICAL TECHNOLOGY CO., LTD
4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue, High-Tech Development Zone
23008 Hefei, China

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SPECIFICITY: 99.8%(95%CI: 94.4%-99.9%)


PRODUCT FEATURES

- ◆ Room temperature storage.
- ◆ No need instrument, get results in 15-20 minutes.
- ◆ Identify acute or early infection
- ◆ No reduction in sensitivity test against the UK, South African and Brazilian variant.


TEST PROCEDURE

1. Specimen Collection

2. Specimen Preparation

3. Testing


Hold the extraction tube vertically and add two drops of the test specimens into the specimen well (S). Start the timer. Interpret the results within 15-20 minutes, and the results after 30 minutes are no longer valid.

4. Interpretation of test results

15-20 minutes

Positive Negative Invalid Invalid

ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD.

【Address】 4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue, High-Tech Development Zone, Hefei 230088, Anhui, China
【Zip code】 230088

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ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD
 4th Floor,D-1# Zone, Pearl Industrial Park, 106 Innovation
 Avenue,High-Tech Development Zone,230088 Hefei,

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1. Introduction

DceepBlue Medical is a national high-tech enterprise specializing in the research and development, production and sales of POCT rapid diagnosis. Its main products include infectious diseases, drugs, gynecology, pregnancy and childbirth testing and other categories. It plans to build the world's first deep blue POCT testing brand. The company has an advanced rapid diagnostic reagent production line, applying the latest clinical molecular diagnostic technology, the product has the advantages of simple operation, accurate and reliable, timely and convenient, and its technical performance has reached the international advanced level. The company has passed the quality management system audit of the Food and Drug Administration and the European Union, and established the quality management system according to the YY/T0287-2017 idt ISO13485:2016 standard, and obtained the European Union EN ISO 13485:2016 quality system certification and CE product certification.

DeepBlue Medical actively carries out cooperation in production, education, and research. It has established extensive technical exchanges and cooperation with universities in the United States and Canada, and internally with the University of Science and Technology of China, Hefei University of Technology, China Pharmaceutical University, Wenzhou University, and various diseases The control center, etc. establish a technology platform and cooperation center to quickly transform scientific research results into productivity. One of the founders has a deep government background and contacts in the pharmaceutical industry, and has studied and flexibly used China's medical policies; the second founder has strong corporate management and operational capabilities and levels in the pharmaceutical industry, and has served as an executive in a number of listed companies. The third founder has a strong technical background, returned scholars, doctoral supervisors, and has a very high academic level in the fields of genetic engineering, molecular diagnosis, rapid diagnosis, disease diagnosis, tumor detection, and food safety. Won the first prize of "Science and Technology Inspection Award" twice by the General Administration of Quality Supervision, Inspection and Quarantine (2010, 2013); published heavyweight papers in internationally renowned journals for many times, and has a number of invention patents.

DeepBlue Medical Device Production License

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CERTIFICATE

No. Q5 18 03 03706 001

Holder of Certificate: ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD.

4th Floor, D-1# Zone
Pearl Industrial Park
106 Innovation Avenue, High-Tech Development Zone
230088 Hefei, Anhui
PEOPLE'S REPUBLIC OF CHINA

Facility(ies):

ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD.
4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue, High-Tech Development Zone, 230088 Hefei, Anhui, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: Design and Development, Production and Distribution of In Vitro Diagnostic Reagents by Colloidal Gold and Enzyme Chemical Reaction Method, Medical Ultrasonic Couplant, Acetowhite Solution, Epithelial Tissue Staining Solution, Rapid Test for Vaginitis (Polyamines) and Cell Preservation Solution

Applied Standard(s):

EN ISO 13485:2016
Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

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DeepBlue Export Product Registration List

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List of antigen tests for direct pathogen detection of coronavirus SARS-CoV-2

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Exemption for COVID 19 Self test in Czech Republic 3

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Medicines & Healthcare products
Regulatory Agency

Medicines & Healthcare products
Regulatory Agency

10 South Colonnade
Canary Wharf
London
E14 4PU
United Kingdom

+44 (0) 20 3080 6000
gov.uk/mhra

lotus global co ltd
23
Maine Street
Reading
Reading
RG2 6AG
England, United Kingdom

26 February 2021

Dear Xueli WEI

We are pleased to confirm that the application to register or update an existing registration for the following manufacturer, which you submitted on **26 February 2021** has been reviewed:

Application reference: **2021022601194338**

Manufacturer organisation: **ANHUI DEEPBLUE MEDICAL TECHNOLOGY CO.,LTD**

Address:

**4th Floor, D-1# Zone, Pearl Industrial Park, 106 Innovation Avenue, High-Tech Development Zone
Anhui
Hefei
230088
China**

Manufacturer registration status: **Registered**

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| GMDN term | Status | MHRA comment |
|--|------------|--------------|
| 64787 - SARS-CoV-2 antigen IVD, kit, immunochromatographic test (ICT), rapid | Registered | |

Please note this letter **does not** represent any form of accreditation, certification or approval by the UK Competent Authority.

If you stop placing devices on the market or if you are not complying with the Regulations, you should inform us so that we can amend our records. You should be aware that it is an offence to place on the market UKCA or CE marked devices that do not comply with the regulations.

Please inform us of the following chargeable changes:

- 1. company/organisation information e.g. name and address**
- 2. additional devices (GMDN code or term)**

Please also use the Devices Online Registration Database (DORS) to tell us of the following changes e.g. removal/discontinuation of a device (GMDN) or product from your registration record, change of contact person, telephone number and/or email address, for which payment of our statutory fee does not apply.

Please note that the name and address of manufacturer, UK Responsible Person or Authorised Representative (Northern Ireland only) and devices that have been registered will be published on our [Public Access Registration Database](#) (PAR). This applies to non-in vitro diagnostic devices only.

The account number for your company/organisation is **0000012384**.

Yours sincerely,

Ngozi Onyeukwu
Device registrations service
Devices division
MHRA

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Anhui Deepblue Medical Technology Co.,Ltd.

Website:<http://www.dbluemedical.com/>

Address:4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue,High-Tech Development Zone 230088 Hefei,Anhui,China

COVID 19 TESTING – LFT STATUS

8 January 2021
Fengling Chen
Anhui Deepblue Medical Technology Co. Ltd

TC/HF/20

Dear Fengling

In light of the public holidays around Christmas and New Year and associated short pause in validation activities during this period, we thought it would be useful to provide a summary of the validation level which your LFDs have reached.

| Name of Product | Sample Method | Validation Level passed |
|--|---------------------------------|---|
| Anhui Deepblue COVID-19 (SARS-CoV-2) Antigen Test Kit (Colloidal Gold) | Nasopharyngeal or Oropharyngeal | All phases up to and including Phase 3B |

CONFIDENTIALITY

We remind you of your ongoing obligations of confidentiality as set out in the non-disclosure agreement recently signed by you dated 8 January 2021.

We thank you for your continued interest in supporting our COVID-19 response.

Yours sincerely,

Oxford University and PHE

<https://www.medrxiv.org/content/10.1101/2021.01.13.21249563v1.full>

<https://www.medrxiv.org/content/10.1101/2021.01.13.21249563v1.full>

medRxiv THE PREPRINT SERVER FOR HEALTH SCIENCES

CSH Cold Spring Harbor Laboratory BMJ Yale

HOME | ABOUT | SUBMIT | NEWS & NOTES | ALERTS / RSS

Search

Advanced Search

[View current version of this article](#)

[Comment on this paper](#)

COVID-19: Rapid Antigen detection for SARS-CoV-2 by lateral flow assay: a national systematic evaluation for mass-testing

Tim Peto, UK COVID-19 Lateral Flow Oversight Team
doi: <https://doi.org/10.1101/2021.01.13.21249563>

This article is a preprint and has not been peer-reviewed [what does this mean?]. It reports new medical research that has yet to be evaluated and so should not be used to guide clinical practice.

Abstract Full Text Info/History Metrics Preview PDF

Abstract

Background New lateral flow device (LFD) viral antigen immunoassays have been developed by commercial and research organisations around the world as diagnostic tests for SARS-CoV-2 infection. To support decisions by the UK Government on potential scale-up of mass population testing, we have at their request evaluated the diagnostic performance of a significant number of point-of-care rapid SARS-CoV-2 LFDs.

Methods 132 LFDs were initially reviewed by a Department of Health and Social Care team, part of the UK government, from which 64 were selected for further evaluation. Standardised laboratory evaluations, and for those that met the published criteria, field testing in the Falcon-C19 research study and UK pilots were performed (UK COVID-19 testing centres, hospital, schools, armed forces).

Results 4/64 LFDs so far have desirable performance characteristics from independent laboratory studies and early preliminary field evaluations (Orient Gene, Deepblue and Innova SARS-CoV-2 Antigen Rapid Qualitative Test), of which one underwent extended clinical assessment in field studies (Innova). 8951 Innova LFD

Posted January 15, 2021.

[Download PDF](#) [Email](#)
[Author Declarations](#) [Share](#)
[Supplementary Material](#) [Citation Tools](#)
[Data/Code](#)

COVID-19 SARS-CoV-2 preprints from medRxiv and bioRxiv

Subject Area
Infectious Diseases (except HIV/AIDS)

Subject Areas

All Articles

- Addiction Medicine
- Allergy and Immunology
- Anesthesia
- Cardiovascular Medicine
- Dentistry and Oral Medicine
- Dermatology
- Emergency Medicine
- Endocrinology (including Diabetes Mellitus and Metabolic Disease)
- Epidemiology

Anhui Deepblue Medical Technology Co.,Ltd.

Website: <http://www.dbluemedical.com/>

Address: 4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue, High-Tech Development Zone 230088 Hefei, Anhui, China

Phase 2

As part of Phase 2 evaluations, 9,692 LFD tests were performed at PHE Porton Down across the 64 candidate devices as of the 3rd December 2020. 5 LFDs had a kit failure rate above the pre-specified threshold for exclusion (>10%), 17 kits had a false-positive rate below the pre-defined specificity threshold (<97%) and 28 kits a false-negative rate below the LOD threshold (<60% at 10² pfu/m). In total, across all three criteria, nineteen kits performed at a level in accordance with the UK Lateral Flow Oversight Group's *a priori* "prioritisation criteria". All nineteen kits also passed cross-reactivity analyses against seasonal human coronaviruses.

Phase 3

To date, eight LFDs have passed Phase 3a evaluation, namely: *Innova SARS-CoV-2 Antigen Rapid Qualitative Test* (Innova), *Zhejiang Orient Gene Biotech Co. Coronavirus Ag Rapid Test Cassette* (Swab) (Orient Gene), *Anhui Deepblue Medical Technology COVID-19 (Sars-CoV-2) Antigen Test kit (Colloidal Gold)* (Deepblue), *Fortress Diagnostics Coronavirus Ag Rapid Test* (Fortress), *Roche SD Biosensor Standard Q COVID-19 Ag Test (SD Bio swab)*, *Surescreen Diagnostics SARS-CoV-2 Antigen Rapid Test Cassette (Nasopharyngeal swab)* (Surescreen) and LFD x (the manufacturer had not given consent to be named). (Supplementary Table 1). Three LFDs did not pass 3a evaluation and the remaining LFDs are currently undergoing evaluation. Four LFDs (Deepblue, Innova, Orientgene, LFD x) have passed Phase 3b evaluation (Table 1, Supp Figure 1), one LFD did not pass and the remainder have not been evaluated.

| Viral Load | Average Ct | Innova Number tested/number positive (%) | LFD x Number tested/number positive (%) | Orient Gene Number tested/number positive (%) | Deepblue Number tested/number positive (%) |
|----------------|------------|--|---|---|--|
| >10million | <18 | 5/5 (100) | 1/1 (100) | - | 3/3 (100) |
| 1-10 million | 18-21.5 | 23/23 (100) | 12/13 (92) | 17/17 (100) | 19/19 (100) |
| 0.1-1 million | 21.5-25 | 52/54 (96) | 19/21 (91) | 18/18 (100) | 43/44 (98) |
| 10,000-100,000 | 25-28 | 37/42 (88) | 13/13 (100) | 18/19 (95) | 38/38 (100) |
| 1,000-10,000 | 28-31 | 25/33 (76) | 17/19 (90) | 14/18 (78) | 18/29 (62) |
| 100-1,000 | 31-34.5 | 11/33 (33) | 10/26 (39) | 11/19 (58) | 8/36 (22) |
| <100 | >34.5 | 2/7 (29) | 1/6 (17) | 0/4 (0) | 0/8 (0) |
| Overall | na | 155/197 (79) | 73/99 (74) | 78/95 (82) | 129/177 (73) |

Table 1. Results of the Phase 3b evaluations showing viral antigen detection/sensitivity of four LFD tests using dry-swab samples from community sampling. Tests were performed by laboratory scientists. Ct – cycle threshold on RT-PCR.

Public Health
England

SARS-CoV-2 Inactivation Testing: Interim Report

| | |
|---|-----------------|
| Report identifier | HCM/CoV2/053/v1 |
| Report date | 30 October 2020 |
| Undertaken by High Containment Microbiology, NIS Laboratories, National Infection Service, Public Health England N.B. This is an interim report and may be updated as further results are obtained | |

| Product/treatment details | |
|---|--|
| Product/treatment | COVID-19 (SARS-CoV-2) Antigen Test Kit Extraction Reagent |
| Manufacturer | Anhui Deepblue Medical Technology Co. |
| Product code | Not known |
| Manufacturer's recommended ratio of sample to product | Swab to be added directly to tube containing approximately 300µl extraction buffer |

| Sample details | |
|--|--|
| Sample type tested | Tissue culture fluid containing 5% (v/v) foetal calf serum |
| Virus strain tested | SARS-CoV-2 England 2 |
| Ratio of spiked virus stock to sample matrix | Not applicable; tissue culture fluid used undiluted |

Report identifier and version number: HCM/CoV2/053/v1

Report date: 30 October 2020

Page 1 of 4

UNCONTROLLED WHEN PRINTED

Anhui Deepblue Medical Technology Co.,Ltd.

Website:<http://www.dbluemedical.com/>

Address:4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue,High-Tech Development Zone 230088 Hefei,Anhui,China

| Experimental conditions | |
|--------------------------------------|---|
| Ratio of sample to product tested | 1 volume sample to 3 volumes product |
| Contact time/s | 1 minute; 5 minutes; 10 minutes |
| Temperature of incubation | Room temperature |
| Brief description of tests performed | Triplicate samples were treated with test buffer for indicated contact time/s or mock-treated in triplicate with an equivalent volume of PBS. Samples were immediately titrated on Vero E6 cells to establish virus titre. This test is quantitative and reports the titre of virus in each treatment condition in TCID ₅₀ per ml. Reduction in virus titre following treatment is given as the difference between the mean log ₁₀ TCID ₅₀ /ml for treated conditions and the PBS control. |

| Table of results | | |
|----------------------------------|--|---|
| | Mean virus titre in log ₁₀ TCID ₅₀ /ml [95% confidence interval] | Titre reduction in log ₁₀ TCID ₅₀ /ml [95% confidence interval] |
| PBS-treated | 6.1 [5.8-6.3] | - |
| Test buffer-treated (1 minute) | 6.1 [5.8-6.4] | 0.0 [-0.4-0.3] |
| Test buffer-treated (5 minutes) | 6.2 [5.9-6.6] | 0.0 [-0.6-0.3] |
| Test buffer-treated (10 minutes) | 6.0 [5.7-6.3] | 0.1 [-0.3-0.5] |

Report identifier and version number: HCM/CoV2/053/v1

Report date: 30 October 2020

Page 2 of 4

UNCONTROLLED WHEN PRINTED

Interpretation

Treatment with Deepblue COVID-19 Antigen Test Kit Extraction Reagent had no effect on SARS-CoV-2 titre following a 1, 5 or 10 minute treatment time.

This buffer should not be relied upon to inactivate infectious samples.

This test has been performed using tissue culture fluid. The effectiveness of this treatment against SARS-CoV-2 may vary when used to inactivate clinical samples or other types of sample matrix. Any results of inactivation testing using other sample matrices will be released as they become available.

Inactivation reagents should not be assumed to be 100% effective against SARS-CoV-2.

Suitability of products and treatments for inactivation of other pathogens has not been evaluated in this study.

All COVID-19 laboratory testing workflows must be subjected to suitable and sufficient risk assessment, with consideration given to any inactivation step. Risk assessments should be reviewed regularly as new information on the inactivation of SARS-CoV-2 becomes available.

The impact of chosen inactivation method on the sensitivity of subsequent SARS-CoV-2 detection should also be assessed locally.

Anhui Deepblue Medical Technology Co.,Ltd.

Website:<http://www.dbluemedical.com/>

Address:4th Floor D-1# Zone, Pearl Industrial Park 106 Innovation Avenue,High-Tech Development Zone 230088 Hefei,Anhui,China

Disclaimer

PHE's evaluations of commercial products and treatments for inactivating SARS-CoV-2 have been carried out primarily for PHE's own internal use and the reports of such evaluations are shared solely for readers information; PHE does not in any way recommend any particular product for virus inactivation; and PHE shall not be responsible for the choice of product or treatment for virus inactivation, and it is the responsibility of the testing laboratory to ensure that any such product or treatment implemented has undergone the necessary verification and validation; and PHE shall not be liable, to the greatest extent possible under any applicable law, for any claim, loss or damage arising out of or connected with use of this and related reports and choice of virus inactivation products or treatments.

PHE is an Executive Agency of the Department of Health and Social Care. Unauthorised use of the PHE name and/or logo is prohibited.

Summary of revisions

Version 1: New document

Queries regarding this report or HCM inactivation testing should be directed to
HCMgroup@phe.gov.uk

Report identifier and version number: HCM/CoV2/053/v1

Report date: 30 October 2020

Page 4 of 4

UNCONTROLLED WHEN PRINTED