



SARS-CoV-2 Rapid Antigen Test

COLLOIDAL GOLD METHOD

The SARS-CoV-2 Rapid Antigen Test (Colloidal Gold Method) is intended for the qualitative detection of the SARS-CoV-2 nucleocapsid protein in human nasopharyngeal or nasal swab sampled from individuals suspected of COVID-19.

- ✓ **Simple**
- ✓ **Rapid**
- ✓ **Convenient**



PRODUCT FEATURES



Sensitivity



Specificity



Test time

METHOD OF USE

Nozzle Cap with Protective Cover



Step 1

Squeeze



Step 2

5 times



Step 3

1 minute




Step 4

SARS-CoV-2 Antigen Test Cassette



Step 5


RESULT INTERPRETATION



POSITIVE



NEGATIVE



INVALID

KIT COMPONENTS



SARS-CoV-2 Antigen Test Cassette



Lysis Buffer



Nozzle Cap with Protective Cover



Disposable Sampler

1 test / pouch × 25 pouches

0.6 ml / tube × 25 tubes

25 devices / bag × 1 bag

25 devices/bag × 1 bag



SCAN ME